Package 'httrlib'

September 19, 2023

Title pipeline for automated in vitro screening platform

Version 0.7.3

Description US EPA has designed a rapid and automated in vitro screening platform using the TempO-seq targeted RNA-seq assay to profile chemical bioactivity across a range of concentrations and cell types. The scale and complexity of these data (measurements of ~20,000 genes across thousands of samples) has necessitated the development of a bioinformatics pipeline to ensure standardized and reproducible results across these screening efforts. This package includes R and Python code for rapid and robust processing of TempO-seq data, leveraging existing open-source tools wherever possible (e.g. HISAT2, samtools, DE-Seq2). This pipeline also includes a battery of QC metrics developed specifically for TempO-seq data generated in large-scale automated experiments, and additional functions to store and access multiple levels of analysis results in an intuitive NoSQL database (MongoDB). The package is intended for both internal EPA/ORD use and to allow external researchers to reproduce our results, process their own data, and ultimately facilitate the development of community standards for high-throughput transcriptomic analysis in a regulatory context.

```
License `use_mit_license()`, `use_gpl3_license()` or friends to pick a
     license
Encoding UTF-8
Roxygen list(markdown = TRUE)
RoxygenNote 7.2.3
exports DESeq2_cut
     addQCplotCols
      as.mongo.date
     countGini
     countQC
     countQCflag
     countStats
      degFrame
      deleteByID
      filterProbesByMean
      filterSamplesByTotal
      findByID
      findIDs
      flattenDEG
      getAnlName
      getChemTrts
      getCredentials
      getDB
```

2 R topics documented:

getDEGs													
getDocIDs													
getFCmatrix													
getProbeManifest													
getQCdefault													
getTrts													
getWellCounts													
getWellInfo													
insertByID													
insertManyDEG													
insertOneDEG													
iterate													
mongoQuery													
mongoURL													
openMongo													
outerFences													
plotQCjitter													
plotQCviolin													
qcBatch													
runDESeq2													
runDESeq2ForChemC	Cond												
runDESeq2Single													
sameCtrlWells													
sampleID_wrapper													
splitIDs													
trtGrpWells													
tss_deg_count													
tss_12_norm													
update_attenuation_fa	ctors												
updateqcFlags													
validateProbeManifest	t												
validate_httr_well_trt_	schema												
Suggests BiocParallel,													
DESeq2,													
ggplot2, reldist,													
rlist,													
stats, testthat													
Imports data.table,													
jsonlite,													
mongolite,													
foreach													
R tonics document	·ha												
R topics document	eu.												
.onLoad													4
addQCplotCols .													5
allBspLysateGroup													6
allBspRNAGroups													6
allBulkI vesteGrou		• •	 	 	• •	 •	• •	 •	• •	 •	 •	• •	7

allRefDoseGroups
allRefGroups
allRefRNAGroups
allTestGroups
as.mongo.date
bspLysateGroup
bspRNAGroup
bulkLysateGroup
check_collection_remapped
chemTreatGroup
countGini
countQC
countQCflag
countStats
degFrame
deleteByID
DESeq2_cut
extract_and_run
fill_dose_for_well_trt
fill_dose_for_well_trt_wrapper
filterProbesByMean
filterSamplesByTotal
Find
findByID
findIDs
flattenDEG
format_for_db_insert
getAnlName
getChemTrts
getCredentials
getDB
getDEGs
getDocIDs
getFCmatrix
getProbeManifest
getQCdefault
getTrts
getWellCounts
getWellInfo
get_global_option
insertByID
insertManyDEG
insertOneDEG
insert_into_httr_study
mongoURL
openMongo
outerFences
plotQCjitter
plotQCviolin
process_chemInfo

4 .onLoad

	process_chem_columns	49
	qcBatch	49
	readWellTrtFile	50
	refChemGroup	51
	refRNAGroup	52
	runDESeq2	53
	runDESeq2ForChemCond	54
	runDESeq2Single	56
	sameCtrlWells	57
	sampleID_scan_and_update	58
	sampleID_wrapper	59
	splitIDs	60
	tests	60
	treatGroupFromSamples	61
	treatPairGroup	62
	trtGrpWells	64
	tss_deg_count	64
	tss_12_norm	65
	unit_subtest_finish	66
	unit_subtest_ident	66
	unit_subtest_length	67
	unit_subtest_set	68
	unit_subtest_type	68
	unit_subtest_values	69
	unit_subtest_values_1	69
	unit_test	70
	updateqcFlags	70
	update_attenuation_factors	71
	validateProbeManifest	72
Index		73

.onLoad

options set upon package loading Options used for the QC flag callsto skip any of these remove them from httrQCflags httrMinMapdFrac
= Flag wells with mapd_frac < this value, default 0.5 httrMinMapdN
= Flag wells with n_reads_mapd < this value, default 300k httrMinNsig80 = Flag wells with n_sig80 < this value, default 1000 httrMinNcov5 = Flag wells with n_cov5 < this value, default 5000 httrMaxTop10Prop = Flag wells with top10_prop > this value, default 0.1
(but NOT used) httrMaxGini = Flag wells with gini_coef > this value,
default 0.95 httrQCflags = Vector of QC flags, order defines flag priority (first flag that hits will be used) - default is to skip HIGH_TOP10
flag min_colsum = floor value under which sample_id are removed
mean_cnt (integer) = floor value to exclude samples with colsum below this value

addQCplotCols 5

Description

options set upon package loading Options used for the QC flag calls - to skip any of these remove them from httrQCflags httrMinMapdFrac = Flag wells with mapd_frac < this value, default 0.5 httrMinMapdN = Flag wells with n_reads_mapd < this value, default 300k httrMinNsig80 = Flag wells with n_sig80 < this value, default 1000 httrMinNcov5 = Flag wells with n_cov5 < this value, default 5000 httrMaxTop10Prop = Flag wells with top10_prop > this value, default 0.1 (but NOT used) httrMaxGini = Flag wells with gini_coef > this value, default 0.95 httrQCflags = Vector of QC flags, order defines flag priority (first flag that hits will be used) - default is to skip HIGH_TOP10 flag min_colsum = floor value under which sample_id are removed mean_cnt (integer) = floor value to exclude samples with colsum below this value

Usage

```
.onLoad(lib, pkg)
```

addQCplotCols

addQCplotCols Function to adjust well treatment data in a way that is better for QC plotting Reclassifies stype/rna_src column to "scat" column, and converts to a factor and creates log10 scale columns from n_reads and n_reads_mapd, and bad_probe_prop if necessary cols present NOTE: Requires data.table package

Description

addQCplotCols Function to adjust well treatment data in a way that is better for QC plotting Reclassifies stype/rna_src column to "scat" column, and converts to a factor and creates log10 scale columns from n_reads and n_reads_mapd, and bad_probe_prop if necessary cols present NOTE: Requires data.table package

Usage

```
addQCplotCols(data)
```

Arguments

data

(*data.frame* or *data.table*) = Table containing all well treatment and QC info, e.g. as returned by getWellInfo()

Value

```
(data.table) = data with additional columns
```

6 allBspRNAGroups

allBspLysateGroups allBspLysateGroups Construct documents for httr_trt_grp_cmp collection corresponding to each BioSpyder Reference QC Lysate sample comparison (single conc treatment). Loop over all plate groups and create a httr_trt_grp_cmp document for each pair of reference RNAs.

Description

allBspLysateGroups Construct documents for httr_trt_grp_cmp collection corresponding to each BioSpyder Reference QC Lysate sample comparison (single conc treatment). Loop over all plate groups and create a httr_trt_grp_cmp document for each pair of reference RNAs.

Usage

```
allBspLysateGroups(
   DB = getOption("DB_NAME"),
   well = getOption("HTTR_WELL_NAME"),
   ...
)
```

Arguments

```
DB (string) = database name with httr_well and httr_trt_grp_cmp collections

well (string) = Name of collection with individual well treatment data, default is

"httr_well"

... = All additional args passed to chemTreatGroup
```

Value

(list) = List of documents that were inserted into trt_grp collection

```
allBspRNAGroups Construct documents for httr_trt_grp_cmp collection corresponding to each BioSpyder Reference RNA QC sample comparison (single conc treatment). Loop over all plate groups and create a httr_trt_grp_cmp document for each pair of reference RNAs.
```

Description

allBspRNAGroups Construct documents for httr_trt_grp_cmp collection corresponding to each BioSpyder Reference RNA QC sample comparison (single conc treatment). Loop over all plate groups and create a httr_trt_grp_cmp document for each pair of reference RNAs.

```
allBspRNAGroups(
   DB = getOption("DB_NAME"),
   well = getOption("HTTR_WELL_NAME"),
   ...
)
```

allBulkLysateGroups 7

Arguments

```
DB (string) = database name with httr_well and httr_trt_grp_cmp collections

well (string) = Name of collection with individual well treatment data, default is

"httr_well"

= All additional args passed to chemTreatGroup
```

Value

(list) = List of documents that were inserted into trt_grp collection

```
allBulkLysateGroups
```

allBulkLysateGroups Construct documents for httr_trt_grp_cmp collection corresponding to each bulk lysate QC sample comparison (single conc treatment).

Description

Loop over all plate groups, chem_id (bulk lysate chemicals only) and create a httr_trt_grp_cmp document for each one.

Usage

```
allBulkLysateGroups(
   DB,
   trt_type = "QC sample",
   rna_src = "Bulk Lysate",
   ctrl_chem = "DMSO",
   well = "httr_well",
   ...
)
```

Arguments

```
(string) = name of database with httr_well and httr_trt_grp_cmp collections collections collections

trt_type (str) = Filter stype to get list of plate groups, chem_id

rna_src (str) = Filter rna_src to get list of plate groups, chem_id

ctrl_chem (str) = Skip this chemical when looping over chem_id

well (str) = Name of collection with individual well treatment data, default is "httr_well ... = All additional args passed to chemTreatGroup
```

Details

Parameters:

Value

(list) = List of documents that were inserted into trt_grp collection Store all documents in a list, also track the number that were skipped due to insufficient replicates with passing QC flags:

8 allRefDoseGroups

allRefDoseGroups

allRefDoseGroups Construct documents for httr_trt_grp_cmp collection corresponding to each reference chemical (single or multi conc).

Description

Loop over all plate groups, media, timeh, chem_id, dose_level (reference chemicals only) and create a httr_trt_grp_cmp document for each one. Alternate version of allRefGroups that also loops over and handles reference chemicals with multiple dose levels Note that to use this version, all reference chemical wells in httr_well_trt MUST have dose_level defined (even single conc reference chemicals) For studies with only single conc reference chemicals that do NOT have a dose_level field, use allRefGroups instead (e.g. for the original MCF-7 screens) Also, unlike allRefGroups, this will automatically append plate group ID as part of the trt_name when creating docs

Usage

```
allRefDoseGroups(
   DB,
   trt_type = "reference chemical",
   trt_prop_fields = c("chem_id", "conc", "conc_unit", "dose_level", "stype"),
   grp_id_opts = c("pg"),
   well,
   exp_doses = NULL,
   ...
)
```

Arguments

```
DB
                  (pymongo.MongoClient) = Open connection to database with httr_well and httr_trt_grp_cmp
                  collections
                  (str) = Filter stype to get list of plate groups, media, timeh, and chem_id
trt_type
trt_prop_fields
                  (list of str) = Fields in httr_well collection that should match up for all trt_wells
                  and the singular value should be propagated to field of same name in httr_trt_grp_cmp
                  collection. Default is to propagate chem_id, conc, conc_unit, dose_level, and
                  stype fields
                  (list) = Options to modify trt_grp_id generation, "pg" is usually sufficient for
grp_id_opts
                  reference chemicals to create distinct group IDs for each plate group
                  (str) = Name of collection with individual well treatment data, default is "httr well"
well
exp_doses(list
                  of str) = Expected number of doses for each reference chemical. This will sup-
                  press some debug messages and only warn if the number of doses does not
                  match. For example, c("GEN" = 8, "SIRO" = 8, "TSA" = 1)
                  (httrplcore.PipelineLogger) = Log handler
log
                  = All additional args passed to chemTreatGroup
**kwargs
```

Details

Parameters:

allRefGroups 9

Value

(list of dict) = List of documents that were inserted into trt_grp collection

allRefGroups	allRefGroups Function to loop over all plate groups, media, and times and then call refChemGroup Construct documents for httr_trt_grp_cmp collection corresponding to each reference chemical (single conc).
--------------	--

Description

Loop over all plate groups, media, timeh, chem_id (reference chemicals only) and create a httr_trt_grp_cmp document for each one.

Usage

```
allRefGroups(DB, trt_type = "reference chemical", well = "httr_well", ...)
```

Arguments

DB	(string) = name of database with httr_well and httr_trt_grp_cmp collections collections collections
trt_type	(str) = Filter stype to get list of plate groups, media, timeh, and chem_id
well	(str) = Name of collection with individual well treatment data, default is "httr_well" = All additional args passed to chemTreatGroup

Details

Parameters:

Value

(list) = List of documents that were inserted into trt_grp collection Store all documents in a list, also track the number that were skipped due to insufficient replicates with passing QC flags:

```
allRefRNAGroups Construct documents for httr_trt_grp_cmp collection corresponding to each reference RNA QC sample comparison (single conc treatment).
```

Description

Loop over all plate groups and create a httr_trt_grp_cmp document for each pair of reference RNAs.

```
allRefRNAGroups(DB, well = "httr_well", ...)
```

10 allTestGroups

Arguments

DB	(string) = name of database with httr_well and httr_trt_grp_cmp collections col-
----	--

lections collections

well (str) = Name of collection with individual well treatment data, default is "httr_well"

... = All additional args passed to chemTreatGroup

Details

Parameters:

Value

(list) = List of documents that were inserted into trt_grp collection Store all documents in a list, also track the number that were skipped due to insufficient replicates with passing QC flags:

allTestGroups a

allTestGroups

Description

Construct documents for httr_trt_grp_cmp collection corresponding to each test chemical and dose level treatment.

Usage

```
allTestGroups(DB, well = "httr_well", exp_doses, ...)
```

Arguments

DB (string) = name of database with httr_well and httr_trt_grp_cmp collections col-

lections collections

well (str) = Name of collection with individual well treatment data, default is "httr_well"

exp_doses (int) = Expected number of doses for each test chemical - this will suppress some

debug messages and only warn if the number of doses does not match ... = All

additional args passed to chemTreatGroup

Details

Loop over all plate groups, media, timeh, chem_id, and dose_level and create a httr_trt_grp_cmp document for each one.

Parameters:

Value

(list) = List of documents that were inserted into trt_grp collection Store all documents in a list, also track the number that were skipped due to insufficient replicates with passing QC flags:

as.mongo.date 11

as.mongo.date Convert "POSIXct" to a list object with one member: $`$date'=(ISODate\ str)$

Description

Given an R datetime object, e.g. as returned by Sys.time(), convert to a structure appropriate for conversion to JSON and mongo insert

Usage

```
as.mongo.date(x)
```

Arguments

```
x (POSIXct) = R datetime object, e.g. as returned by Sys.time()
```

Value

(list) = Has one member: '\$date'=(ISODate str), converting this to JSON with auto_unbox=T will insert this as a proper date

bspLysateGroup

bspLysateGroup Construct a document for httr_trt_grp_cmp collection corresponding to the BioSpyder reference QC lysate sample comparison (A vsB). Given a specific plate group, plus any optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This function is just a convenenience wrapper to chemTreatGroup with modified defaults. Note that unlike the original MCF7 bulk lysates, these samples are annotated more like the reference RNA samples.

Description

bspLysateGroup Construct a document for httr_trt_grp_cmp collection corresponding to the BioSpyder reference QC lysate sample comparison (A vsB). Given a specific plate group, plus any optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This function is just a convenenience wrapper to chemTreatGroup with modified defaults. Note that unlike the original MCF7 bulk lysates, these samples are annotated more like the reference RNA samples.

```
bspLysateGroup(
  DB = getOption("DB_NAME"),
  pg_id,
  trt_name = "BSP_LYSATE_B",
  ctrl_name = "BSP_LYSATE_A",
  trt_type = "QC sample",
  ctrl_type = "QC sample",
  trt_src = "Bulk Lysate",
```

12 bspRNAGroup

```
ctrl_src = "Bulk Lysate",
  trt_prop_fields = c(),
  both_prop_fields = c("stype", "rna_src", "pg_id", "block_id"),
  ctrl_desc_field = "trt_name",
  grp_id_opts = c("vs", "pg"),
  ...
)
```

Arguments

DB (string) database name with httr_well and httr_trt_grp_cmp collections (string) = Filter both trt and ctrl wells, match to pg_id in httr_well, required for pg_id ref RNA groups (string) = Match to stype in httr_well, defaults to "QC sample" trt_type ctrl_type (string) = Match to stype in httr_well, defaults to "QC sample" (string) = Match to rna_src in httr_well, defaults to "BSP_RNA_B" trt_src (string) = Match to rna_src in httr_well, defaults to "BSP_RNA_A" ctrl_src trt_prop_fields (list) = Fields in httr well collection that should match up for all trt wells and the singular value should be propagated to field of same name in httr_trt_grp_cmp collection. Default is to propagate just the stype field (exclude all chemical and dose information). both_prop_fields (list) = Fields in httr_well collection that should match up for all trt_wells AND ctrl_wells, and the singular value should be propagated to to field of same name in httr_trt_grp_cmp collection. Default is to propagate pg_id and block_id only for ref RNA comparisons ctrl_desc_field (string) = Field to propagate from ctrl_wells to ctrl field in httr_trt_grp_cmp document, defaults to "trt_name" for ref RNA comparisons (list) = Options to modify trt_grp_id generation, "vs" will capture the pairwise grp_id_opts comparison type, and "pg" is usually sufficient for ref RNA to create distinct group IDs for each plate group

All additional args passed to treatGroupFromSamples

Value

(list) = Document that was inserted into trt_grp collection

bspRNAGroup

bspRNAGroup Construct a document for httr_trt_grp_cmp collection corresponding to the BioSpyder purified reference RNA QC sample comparison (A vsB). Given a specific plate group, plus any optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This function is just a convenenience wrapper to chemTreatGroup with modified defaults.

bspRNAGroup 13

Description

bspRNAGroup Construct a document for httr_trt_grp_cmp collection corresponding to the BioSpyder purified reference RNA QC sample comparison (A vsB).

Given a specific plate group, plus any optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This function is just a convenenience wrapper to chemTreatGroup with modified defaults.

Usage

```
bspRNAGroup(
  DB = getOption("DB_NAME"),
  pg_id,
  trt_name = "BSP_RNA_B",
  ctrl_name = "BSP_RNA_A",
  trt_type = "QC sample",
  ctrl_type = "QC sample",
  trt_src = "Reference RNA",
  ctrl_src = "Reference RNA",
  trt_prop_fields = c(),
  both_prop_fields = c("stype", "rna_src", "pg_id", "block_id"),
  ctrl_desc_field = "trt_name",
  grp_id_opts = c("vs", "pg"),
  ...
)
```

Arguments

```
DB (db name) = Open connection to database with httr_well and httr_trt_grp_cmp collections

pg_id (string) = Filter both trt and ctrl wells, match to pg_id in httr_well, required for ref RNA groups

trt_type (string) = Match to stype in httr_well, defaults to "QC sample"

ctrl_type (string) = Match to stype in httr_well, defaults to "QC sample"

trt_src (string) = Match to rna_src in httr_well, defaults to "BSP_RNA_B"

ctrl_src (string) = Match to rna_src in httr_well, defaults to "BSP_RNA_A"

trt_prop_fields

(vector) = Fields in httr_well collection that should match up for all trt_wells and
```

(vector) = Fields in httr_well collection that should match up for all trt_wells and the singular value should be propagated to field of same name in httr_trt_grp_cmp collection. Default is to propagate just the stype field (exclude all chemical and dose information).

both_prop_fields

(vector) = Fields in httr_well collection that should match up for all trt_wells AND ctrl_wells, and the singular value should be propagated to to field of same name in#' @param httr_trt_grp_cmp collection. Default is to propagate pg_id and block_id only for ref RNA comparisons

ctrl_desc_field

(string) = Field to propagate from ctrl_wells to ctrl field in httr_trt_grp_cmp document, defaults to "trt_name" for ref RNA comparisons

grp_id_opts (vector) = Options to modify trt_grp_id generation, "vs" will capture the pairwise comparison type, and "pg" is usually sufficient for ref RNA to create distinct group IDs for each plate group 14 bulkLysateGroup

Value

(list) = Document that was inserted into trt_grp collection

bulkLysateGroup

bulkLysateGroup Construct a document for httr_trt_grp_cmp collection corresponding to a specific bulk lysate chemical treatment.

Description

Given a specific plate group, plus any optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This function is just a convenenience wrapper to chemTreatGroup with modified defaults.

Usage

```
bulkLysateGroup(
   DB,
   pg_id,
   trt_chem = "TSA",
   trt_dose = NULL,
   trt_type = "QC sample",
   ctrl_type = "QC sample",
   ctrl_chem = "DMSO",
   rna_src = "Bulk Lysate",
   trt_prop_fields = c("chem_id", "conc", "conc_unit", "stype"),
   both_prop_fields = c("pg_id", "block_id"),
   ctrl_desc_field = "chem_id",
   grp_id_opts = c("pg"),
   ...
)
```

Arguments

```
DB
                  (string) = name of database with httr_well and httr_trt_grp_cmp collections col-
                  lections collections
                   (str) = Filter both trt and ctrl wells, match to pg_id in httr_well, required for bulk
pg_id
                  lysate groups
                  (str) = Match to chem id in httr well, defaults to "TSA"
trt chem
                  (int) = Match to dose_level in httr_well - Excluded by default
trt_dose
                  (str) = Match to stype in httr_well, defaults to "QC sample"
trt_type
                  (str) = Match to stype in httr_well, defaults to "QC sample"
ctrl_type
ctrl_chem
                  (str) = Match to chem_id in httr_well, defaults to "DMSO"
rna_src
                  (str) = Match to rna src in httr well for both trt and ctrl wells, defaults to "Bulk
                  Lysate"
trt_prop_fields
```

(vector) = Fields in httr_well collection that should match up for all trt_wells and the singular value should be propagated to field of same name in httr_trt_grp_cmp collection. Default is to propagate chem_id, conc, conc_unit, and stype fields (exclude dose_level here).

both_prop_fields

(vector) = Fields in httr_well collection that should match up for all trt_wells AND ctrl_wells, and the singular value should be propagated to to field of same name in httr_trt_grp_cmp collection. Default is to propagate pg_id and block_id only for bulk lysate comparisons

ctrl_desc_field

(str) = Field to propagate from ctrl_wells to ctrl field in httr_trt_grp_cmp document, defaults to "chem_id" for bulk lysate comparisons

grp_id_opts

(vector) = Options to modify trt_grp_id generation, "pg" is usually sufficient for bulk lysate to create distinct group IDs for each plate group ... = All additional args passed to treatGroupFromSamples

Details

Parameters:

Value

(list) = Document that was inserted into trt_grp collection Pass to chemTreatGroup:

```
check_collection_remapped
```

check_collection_remapped Check if any collections are re-mapped, fill in defaults for everything else

Description

check_collection_remapped Check if any collections are re-mapped, fill in defaults for everything else

Usage

```
check_collection_remapped(collections)
```

Arguments

collections (named character vector) = vector to re-map collection names used, where name=default collection; value=new collection name

Value

collections (*named character vector*) = vector to re-map collection names used, where name=default collection; value=new collection name

16 chemTreatGroup

chemTreatGroup

chemTreatGroup Function that generates the trt_grp_cmp for a specific test chemical treatment vs DMSO Construct a document for httr_trt_grp_cmp collection corresponding to a specific chemical and dose level treatment.

Description

Given a specific chemical and dose, plus optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. Note, if filtering reduces either trt or ctrl group to < min_reps samples this treatment group will be skipped.

Usage

```
chemTreatGroup(
  DB,
  trt_chem,
  trt dose,
  trt_type = "test sample",
  ctrl_type = "vehicle control",
  ctrl_chem = "DMSO",
  cell_type = NULL,
  pg_id = NULL_i
  media = NULL,
  timeh = NULL,
  trt_src = NULL,
  ctrl_src = NULL,
  qc_flags = c("OK"),
  min\_reps = 2,
  well = "httr_well",
)
```

Arguments

```
DB
                   (string) = name of database with httr_well and httr_trt_grp_cmp collections col-
                   lections
                   (str) = Match to chem_id in httr_well
trt_chem
trt_dose
                   (int) = Match to dose_level in httr_well (set to None to exclude)
                   (str) = Match to stype in httr_well, defaults to "test sample"
trt_type
ctrl_type
                   (str) = Match to stype in httr_well
                   (str) = Match to chem_id in httr_well
ctrl_chem
pg_id
                   (str) = Filter both trt and ctrl wells, match to pg_id in httr_well, defaults to None
                   (str) = Filter both trt and ctrl wells, match to media in httr_well, defaults to None
media
timeh
                   (str) = Filter both trt and ctrl wells, match to timeh in httr_well, defaults to None
                   (str) = Filter trt wells, match to rna_src field
trt_src
                   (str) = Filter ctrl wells, match to rna_src field
ctrl_src
```

countGini 17

qc_flags	(vector) = Filter both trt and ctrl wells, match to qc_flag field, defaults to OK only
min_reps	(int) = Minimum number of replicates in trt and ctrl groups, respectively - if either is less, will generate an output message and return empty dict
well	$(str) = Name of collection with individual well treatment data, default is "httr_well"$
	= All additional args passed to treatGroupFromSamples

Details

Parameters:

Value

(list) = Document that was inserted into trt_grp collection Build the query for treatment and ctrl wells

countGini

countGini Computing Gini coefficient Use reldist::gini(counts) where counts = probe counts for a single sample This metric DOES depend on whether 0 counts are present or missing from the vector If zero counts are missing, should specify proben to fill in missing 0s

Description

countGini Computing Gini coefficient Use reldist::gini(counts) where counts = probe counts for a single sample This metric DOES depend on whether 0 counts are present or missing from the vector If zero counts are missing, should specify proben to fill in missing 0s

Usage

```
countGini(
  counts,
  proben = getOption("httrProbeN"),
  debug = getOption("debug", default = FALSE)
)
```

Arguments

```
counts (integer) = a vector of counts for a single sample, w/ or w/out 0 counts

proben (integer) = Total number of probes used for Gini coefficients across samples, ignored if NULL

debug (logical) = Whether or not to report extra debug messages
```

Value

(numeric) = Gini coefficient returned from reldist::gini

18 countQC

countQC	countQC Given a document from httr_counts collection, compute all QC stats and flags and format according to httr_counts_qc schema option() function should be used to change defaults for underlying QC functions

Description

countQC Given a document from httr_counts collection, compute all QC stats and flags and format according to httr_counts_qc schema option() function should be used to change defaults for underlying QC functions

Usage

```
countQC(
  counts_doc,
  bad_probes = c(),
  proben = getOption("httrProbeN"),
  calc_gini = getOption("httrStatGini", default = TRUE),
  calc_flag = TRUE,
  min_mapd_frac = get_global_option("min_mapd_frac"),
  min_n_reads_mapd = get_global_option("min_n_reads_mapd"),
  min_n_sig80 = get_global_option("min_n_sig80"),
  min_n_cov5 = get_global_option("min_n_cov5"),
  max_top10_prop = get_global_option("max_top10_prop"),
  max_gini_coef = get_global_option("max_gini_coef"),
  qc_flags = get_global_option("qc_flags"),
  debug = getOption("debug", default = FALSE)
)
```

Arguments

counts_doc	(list) = Follows DB schema for httr_counts, msut include the _id field
bad_probes	(character vector) = List of probe IDs that should be masked
proben	(integer) = Total number of probes on the platform (length bad_probes will be subtracted before passing to gini function)
calc_gini	(<i>logical</i>) = Whether or not to calculate gini, default True but can override if no reldist package installed
calc_flag	$(logical)$ = Whether or not to calculate qc_flag, default True but when set to False all flags will be set to "OK"
debug	(logical) = Whether to report additional debug info

countQCflag 19

countQCflag

countQCflag Given a list structure with all the quantitative fields in the httr_counts_qc schema, determine the qc_flag field Flags are checked in priority order according to qc_flags parameter NOTE: If any required fields are missing from qc_stats, it will lead to an error

Description

countQCflag Given a list structure with all the quantitative fields in the httr_counts_qc schema, determine the qc_flag field Flags are checked in priority order according to qc_flags parameter NOTE: If any required fields are missing from qc_stats, it will lead to an error

Usage

```
countQCflag(
   qc_stats,
   min_mapd_frac = get_global_option("min_mapd_frac"),
   min_n_reads_mapd = get_global_option("min_n_reads_mapd"),
   min_n_sig80 = get_global_option("min_n_sig80"),
   min_n_cov5 = get_global_option("min_n_cov5"),
   max_top10_prop = get_global_option("max_top10_prop"),
   max_gini_coef = get_global_option("max_gini_coef"),
   qc_flags = get_global_option("qc_flags"),
   debug = getOption("debug", default = FALSE)
)
```

Arguments

```
(list) = list of all QC stats computed by countStats and countGini above
qc_stats
min_mapd_frac
                 (numeric) = minimum mapd_frac, flag with LOW_MAPD_FRAC below this
                 cutoff
min_n_reads_mapd
                 (integer) = minimum n reads mapd, flag with LOW READS below this cutoff
                 (integer) = minimum n_sig80, flag with LOW_NSIG80 below this cutoff
min_n_sig80
min_n_cov5
                 (integer) = minimum n_cov5, flag with LOW_NCOV5 below this cutoff
max_top10_prop
                 (numeric) = maximum top10_prop, flag with HIGH_TOP10 above this cutoff
max_gini_coef
                 (numeric) = maximum gini_coef, flag with HIGH_GINI above this cutoff
                 (character vector) = list of flags to apply in priority order
qc_flags
                 (logical) = Whether to report some debug info
debug
```

Value

```
(character) = a single QC flag for the well, "OK" if not thresholds were hit
```

20 degFrame

countStats	countStats Function to compute count-based QC statistics for a single sample
------------	--

Description

countStats Function to compute count-based QC statistics for a single sample

Usage

```
countStats(
  counts,
  ncov = getOption("httrStatNCov", default = 5),
  nsig = getOption("httrStatNsig", default = 0.8),
  topn = getOption("httrStatTopN", default = 10)
)
```

Arguments

counts	(<i>integer</i>) = a vector of counts for a single sample, w/ or w/out 0 counts All of these statistics should be invariant to the probe labels, ordering of the counts, or whether 0 counts are missing
ncov	(<i>integer</i>) = Minimum read count(s) for Ncov statistic (number of probes covered at a minimum read count level)
nsig	(numeric) = Proportion of signal for Nsig statistic (minimum number of probes capturing top X% of signal)
topn	(<i>integer</i>) = Number of probes to use for TopN statistic (proportion of signal captured by top probe)

Value

 $(named\ list)$ = each of the requested statistics, member names now match DB schema

degFrame	degFrame - Take just the degs member of a httr_deg doc and reformat as a data.frame with appropriate classes for each column This a low-level function, is automatically applied by getDEGs when pulling from
	database.

Description

degFrame - Take just the degs member of a httr_deg doc and reformat as a data.frame with appropriate classes for each column This a low-level function, is automatically applied by getDEGs when pulling from database.

```
degFrame(degs)
```

deleteByID 21

Arguments

degs (*list*) - The raw degs member of a doc in httr_deg collection, as returned directly from mongo

Value

(data.frame) - Table with 7 columns matching the underlying deg fields

deleteByID deleteByID

• Delete documents from a collection. Given a list of IDs, delete all matching documents from a collection, or delete ALL documents. The primary purpose is to remove documents from a collection before replacing them when rerun=True in a particular pipeline step.

Description

deleteByID

• Delete documents from a collection. Given a list of IDs, delete all matching documents from a collection, or delete ALL documents. The primary purpose is to remove documents from a collection before replacing them when rerun=True in a particular pipeline step.

Usage

```
deleteByID(
   DB,
   ids = c(),
   id_field = "sample_id",
   delete_all = FALSE,
   debug = getOption("debug", default = FALSE))
```

Arguments

Value

(numeric) = Number of documents that were deleted, will also warn if != length(ids) when delete_all=False

22 extract_and_run

DESeq2_cut	DESeq2_cut- (compartmentalized unit test) - test whether a localled
	generated DESeq2 result matches the the saved result

Description

DESeq2_cut- (compartmentalized unit test) - test whether a localled generated DESeq2 result matches the the saved result

Usage

```
DESeq2_cut(saved_results, treatments)
```

Arguments

```
saved_results

(list) = object with archived DESeq2 results and all necessary objects to recreate that object with local instalation

treatments

(data.frame) = dataframe containing information about the analysis parameters that went into calculating DESeq2 results (shrinkage, normalization)
```

Value

(boolean) = if all tests pass, returns true, if any fail, returns false

```
extract_and_run
```

extract_and_run - Run DESeq2 on a single chemical dose-response series Uses DB access functions to extract all httr_trt_grp_id entries corresponding to a dose-response series, then extracts counts for all wells, filters by min read depth and mean floor, and passes through runDESeq2.

Description

extract_and_run - Run DESeq2 on a single chemical dose-response series Uses DB access functions to extract all httr_trt_grp_id entries corresponding to a dose-response series, then extracts counts for all wells, filters by min read depth and mean floor, and passes through runDESeq2.

```
extract_and_run(
  chem_id,
  db_host,
  db_name,
  pg_id = NULL,
  media = NULL,
  timeh = NULL,
  max_dose_level = 8,
  min_colsum = getOption("min_colsum"),
```

fill_dose_for_well_trt 23

```
mean_cnt = getOption("mean_cnt"),
plate_effect = F,
shrinkage = "normal",
threads = 1,
collections = character(0),
debug = getOption("debug", default = FALSE)
```

Arguments

chem_id	$(character) = ID$ for the chemical to run, matched against chem_id field in httr_trt_grp_cmp
db_host	(character) = Host for DB connection
db_name	(character) = Name of DB
pg_id	(character) = Limit httr_trt_grp_cmp to this plate group, useful if multiple dose response series for same chem_id
media	(character) = Limit httr_trt_grp_cmp to this media type, useful if multiple dose response series for same chem_id
timeh	(integer) = Limit httr_trt_grp_cmp to this timeh, useful if multiple dose response series for same chem_id
max_dose_lev	
	(integer) = Max number of dose levels expected
min_colsum	(integer) = Filter samples by min read total, set to NULL to skip, defaults to $100k$
mean_cnt	(integer) = Filter probes by mean floor, set to NULL to skip, defaults to 5
plate_effect	(logical) = Whether or not to include plate effect in model, defaults to FALSE
shrinkage	(character) = Type of shrinkage to use, defaults to new default, a.k.a. "normal"
threads	(integer) = Number of threads to use for core DESeq2 functions, defaults to 1 (no multi-threading)
collections	(named character vector) = Optional vector to re-map collection names used, where name=default collection; value=new collection name
debug	(logical) = Whether to print debug messages, default: FALSE, overridden by options(debug=)

Value

vector of my_results, chemTrts, trt_grp_name

```
fill_dose_for_well_trt

fill_dose_for_well_trt This function takes a data.table and fills the

dose column with calculated values
```

Description

fill_dose_for_well_trt This function takes a data.table and fills the dose column with calculated values

24 filterProbesByMean

Usage

```
fill_dose_for_well_trt(sampleKey)
```

Arguments

sampleKey:

(data.table) the data.table conforming with the normal welltrt columns

process: for each conc found ordered by pg_id/chem_id assign a corresponding

dose_level

Value

the modified data.table

```
fill_dose_for_well_trt_wrapper
```

fill_dose_for_well_trt_wrapper This function takes a sample_id file and converts it to a data.table and does some checks then calls fill_dose_for_well_trt

Description

fill_dose_for_well_trt_wrapper This function takes a sample_id file and converts it to a data.table and does some checks then calls fill_dose_for_well_trt

Usage

```
fill_dose_for_well_trt_wrapper(sampleID_file)
```

Arguments

sampleKey:

(file) the file conforming with the normal welltrt columns

Value

the dat.frame resulted from fill_dose_for_well_trt

filterProbesByMean filterProbesByMean - Function to prefilter DESeq2 input data by removing probes with low mean count Removes any probes with mean count < mean_cnt

Description

filterProbesByMean - Function to prefilter DESeq2 input data by removing probes with low mean count Removes any probes with mean count < mean_cnt

filterSamplesByTotal 25

Usage

```
filterProbesByMean(
   COUNTS,
   mean_cnt = getOption("mean_cnt"),
   debug = getOption("debug", default = FALSE)
)
```

Arguments

counts (data.frame) = Data frame containing all counts to be analyzed, row.names = probes, colnames = samples

mean_cnt (integer) = Remove any samples with colsum below this value

debug (logical) = Whether to print debug messages, default: FALSE, overridden by

options(debug=...)

Value

(data.frame) = COUNTS after removing the probes with low mean counts

```
filterSamplesByTotal
```

filterSamplesByTotal - Function to prefilter DESeq2 input data by removing low-depth samples Removes any sample with colsum < min_colsum

Description

 $filter Samples By Total - Function \ to \ prefilter \ DESeq2 \ input \ data \ by \ removing \ low-depth \ samples \ Removes \ any \ sample \ with \ colsum < min_colsum$

Usage

```
filterSamplesByTotal(
   COUNTS,
   min_colsum = getOption("min_colsum"),
   debug = getOption("debug", default = FALSE)
)
```

Arguments

```
COUNTS (data.frame) = Data frame containing all counts to be analyzed, row.names = probes, colnames = samples

min_colsum (integer) = Remove any samples with colsum below this value

debug (logical) = Whether to print debug messages, default: FALSE, overridden by options(debug=...)
```

Value

(data.frame) = COUNTS after removing the samples with low colsum

26 findByID

Find

title Find: with DB object as well as fields and potentially query named list, inquire db for matching documents

Description

title Find: with DB object as well as fields and potentially query named list, inquire db for matching documents

Usage

```
Find(DB = NULL, ...)
```

Arguments

DB (mongo object) = Opened connection to specific Mongo Database and Collec-

tion

optional fields: json list of fields to be included in our search other parameter (optional):

query JSON query to use when querying the db

Value

list of documents matching search specified by query/fields

findByID

findByID Find all documents in a collection with optional ID filtering.

Description

Given a vector of IDs (typically either _id or sample_id), return an iterable object (mongo iterator or list) corresponding to all matching documents in a collection. If ids is an empty list, return all documents in the collection. Also has options to limit the result to just the id_field or any other subset of fields, and to dump out the data as a list (stored in local memory) instead of returning the pymongo Cursor object.

```
findByID(
   DB,
   ids = c(),
   id_field = "sample_id",
   id_only = FALSE,
   fields = "{}",
   dump = FALSE,
   debug = getOption("debug", default = FALSE)
)
```

findIDs 27

Arguments

DB	(emphmongo object) = Open connection to specific Mongo Database and Collection
ids	(character vector) = Vector of IDs (usually character) to query on, gets all values if empty, default: empty vector.
id_field	(character) = Which ID field to match against ids, default: sample_id.
id_only	(logical) = If True, return only the values in the ID field being matched against, otherwise return the whole document
fields	(Mongo query, as character) = If id_only = False, this will be used to filter the return fields instead.
dump	(logical) = If True, loop over the initial return object (pymongo.cursor.Cursor) and extract all results to a new list object.
debug	(logical) = If True, dump out some extra debugging messages, can also set using options(debug=True)

Value

mongo iterator, list, or vector) = Either: an iterator object to return results from Mongo (dump=False), a list of all results (dump=True, id_only=False) a vector of IDs (dump=True, id_only=True)

findIDs	findIDs Generic function to get IDs matching any query from any col-
	lection Given an open DB connection or connection parameters, run
	a query and get all distinct values of id_field

Description

findIDs Generic function to get IDs matching any query from any collection Given an open DB connection or connection parameters, run a query and get all distinct values of id_field

Usage

```
findIDs(
   DB = NULL,
   db_host = NULL,
   db_name = NULL,
   collection = NULL,
   id_field = "_id",
   debug = getOption("debug", default = FALSE),
   ...
)
```

Arguments

```
    (mongo object) = If specified, use this open connection to a collection, ignore db_host, db_name, collection
    (character) = If DB is NULL, this specifies the host to open a connection to
    (character) = If DB is NULL, this specifies the DB name to connect to
```

28 format_for_db_insert

collection (character) = If DB is NULL, this specifies the collection to connect to
id_field (character) = What field in the DB to return IDs from, defaults to the generic
mongo_id field

debug (logical) = Whether to print debug messages, default: FALSE, overridden by
options(debug=...) ... = All additional params passed to mongoQuery to form
the query string

Value

(character vector) = The ID values returned by DB\$distinct call

flattenDEG Convert a single httr_deg schema doc to a data.frame This

is a low-level function that is automatically called by getDEGs whenever flatten=T This takes all other fields besides degs and makes them columns in the data.frame anl_opt and update_notes are flattened us-

ing jsonlite::toJSON

Description

flattenDEG Convert a single httr_deg schema doc to a data.frame This is a low-level function that is automatically called by getDEGs whenever flatten=T This takes all other fields besides degs and makes them columns in the data.frame anl_opt and update_notes are flattened using json-lite::toJSON

Usage

```
flattenDEG(deg_doc)
```

Arguments

deg_doc (list) = Single mongo doc conforming to httr_deg schema

Value

(data.frame) = Flatted version of deg_doc

format_for_db_insert

format_for_db_insert description: optionally saves the results provided by prior call to extract_and_run function back to the DB in addition to returning a single data.frame with all results that can be written to a TSV file or saved in Rdata format.

Description

format_for_db_insert description: optionally saves the results provided by prior call to extract_and_run function back to the DB in addition to returning a single data.frame with all results that can be written to a TSV file or saved in Rdata format.

getAnlName 29

Usage

```
format_for_db_insert(
   extract_and_run_return,
   mean_cnt,
   plate_effect,
   shrinkage,
   collections,
   db_host,
   db_name,
   note,
   rerun,
   debug,
   db_insert
)
```

Arguments

extract	and	run	return

vector containing the results computed by prior call to runDESeq2 function and

to be inserted into db

mean_cnt (integer) = Filter probes by mean floor, set to NULL to skip, defaults to 5

 $\verb|plate_effect| (\textit{logical}) = Whether or not to include plate effect in model, defaults to FALSE$

shrinkage (character) = Type of shrinkage to use, defaults to new default, a.k.a. "normal"

collections (named character vector) = Optional vector to re-map collection names used,

where name=default collection; value=new collection name

db_host (character) = Host for DB connection

 db_name (character) = Name of DB

note (*character*) = Overwrite default note fields for insert if this is specified.

rerun (logical) = Whether to re-run the DB insertion if results are in there already

(ideally this should check the DB before running...)

debug (logical) = Whether to print debug messages, default: FALSE, overridden by

options(debug=...)

db_insert (logical) = Whether to insert the results back into the DB, defaults to TRUE

Value

(data.frame) = Contains DESeq2 results for all treatments in the dose-response series

getAnlName	getAnlName Convert analysis options to standardized anl_name for
	httr_deg docs

Description

getAnlName Convert analysis options to standardized anl_name for httr_deg docs

30 getChemTrts

Usage

```
getAnlName(
  mean_cnt = getOption("mean_cnt"),
  plate_effect = F,
  shrinkage = "normal",
   ...
)
```

Arguments

```
mean_cnt (numeric)
plate_effect (character)
shrinkage (character)
```

getChemTrts

getChemTrts

Pull out documents from httr_trt_grp_cmp collection corresponding to a single chemical dose-response series Returns the relevant documents from httr_trt_grp_cmp as a named list, where names come from trt_grp_id field, and are sorted by dose_level

Description

getChemTrts

• Pull out documents from httr_trt_grp_cmp collection corresponding to a single chemical doseresponse series Returns the relevant documents from httr_trt_grp_cmp as a named list, where names come from trt_grp_id field, and are sorted by dose_level

```
getChemTrts(
  DB = NULL,
  chem_id,
  db_host = NULL,
  db_name = NULL,
  collection = "httr_trt_grp_cmp",
  pg_id = NULL,
  media = NULL,
  timeh = NULL,
  max_dose_level = 8,
  debug = getOption("debug", default = FALSE)
)
```

getCredentials 31

Arguments

DB	(mongo object) = If specified, use this open connection to httr_trt_grp_cmp collection, ignore db_host, db_name, collection
chem_id	(character) = Match against chem_id field, this is required
db_host	(character) = If DB is NULL, this specifies the host to open a connection to
db_name	(character) = If DB is NULL, this specifies the DB name to connect to
collection	(character) = If DB is NULL, this specifies the collection to connect to, default: httr_trt_grp_cmp
pg_id	(character) = Optional addl filter on pg_id field, useful if there is more than one series for same chem_id
media	(<i>character</i>) = Optional addl filter on media field, useful if there is more than one series for same chem_id
timeh	(<i>numeric</i>) = Optional addl filter on timeh field, useful if there is more than one series for same chem_id
max_dose_lev	
	(<i>numeric</i>) = Expected max number of doses, will give a warning if more than this number of treatment groups returned, default is 8
debug	(<i>logical</i>) = Whether to print debug messages, default: FALSE, overridden by options(debug=)

Value

(named list) = Names come from trt_grp_id field, each member is a complete httr_trt_grp_cmp doc (list type), members are sorted by dose_level field

lib/db/pc it just qı	ntials This sswords.py I tickly scans wd combo	Rather than	i load enti	re ke	ychains int	o memo	ory,
-------------------------	---	-------------	-------------	-------	-------------	--------	------

Description

getCredentials This function replaces the Keychain class in lib/db/passwords.py Rather than load entire keychains into memory, it just quickly scans ~/.mongopw and ~/.mngdb/passwd for relevant user,passwd combo

Usage

```
getCredentials(host = getOption("httrDefaultHost"), db = NULL)
```

Arguments

```
host (character)
db (mongo object)
```

Value

NULL or a list with members user,passwd (and host,db if specified)

32 getDEGs

getDB

getDB if DB not provided already, opens a connection to mongodb via the mongolite package to the collection provided

Description

getDB if DB not provided already, opens a connection to mongodb via the mongolite package to the collection provided

Usage

```
getDB(DB = NULL, db_host = NULL, db_name = NULL, collection = NULL)
```

Arguments

DB (mongo object) = Opened connection to specific Mongo Database and Collec-

tion or null if not provided

db_host: (character) string representing the db url
db_name: (character) string representing the db name

collection: (character) string representing the collection name we are trying to connect to

Value

(mongo object) a mongo 'environment', i.e. a DB object

getDEGs

getDEGs Pull out DESeq2 results from httr_deg table

Description

getDEGs Pull out DESeq2 results from httr_deg table

```
getDEGs(
  DB = NULL,
  db_host = NULL,
  db_name = NULL,
  collection = "httr_deg",
  single = F,
  flatten = F,
  warn_count = 100,
  warn_stype = T,
  warn_anl = T,
  debug = getOption("debug", default = FALSE),
  ...
)
```

getDocIDs 33

Arguments

DB	(<i>mongo object</i>) = If specified, use this open connection to httr_deg collection, ignore db_host, db_name, collection
db_host	(character) = If DB is NULL, this specifies the host to open a connection to
db_name	(character) = If DB is NULL, this specifies the DB name to connect to
collection	(character) = If DB is NULL, this specifies the collection to connect to, default: httr_deg
single	(logical) = Whether to return a single doc - if TRUE and query matches multiple docs, then only the first match is returned
flatten	(<i>logical</i>) = Whether to flatten all results into a single data frame instead of returning in list structure closer to mongo doc structure
warn_count	(integer) = If query is going to return more than X documents, issue a warning (set to Inf to silence), default: 100
warn_stype	(<i>logical</i>) = When TRUE: If query returns documents with a mix of stype values, issue a warning
warn_anl	(<i>logical</i>) = When TRUE: If query returns documents with a mix of anl_name values, issue a warning
debug	(logical) = Whether to print debug messages, default: FALSE, overridden by options(debug=) = All other params are passed to mongoQuery to build the filtering query

Value

(*list* or *data.frame*) = When single & !flatten, this is a list matching httr_deg schema; when !single & !flatten, this is a list of lists with multiple httr_deg docs; when flatten, this is a data.frame

Description

extract a certain ID field (e.g. sample_id) out of a list of documents This is designed to mimic the function of the same name from lib/db/mongo.py

Usage

```
getDocIDs(docs, id_field = "sample_id")
```

Arguments

docs (list) = List of documents dumped out from MongoDB as a list type id_field (character) = Name of field to pull out of every member of docs

Value

(vector) = Vector of ID values from docs, atomic type depends on the type of ID (usually character)

34 getFCmatrix

getFCmatrix	getFCmatrix - Pull out DESeq2 results for multiple contrasts from
	httr_deg collection and reshape into wide matrix

Description

getFCmatrix - Pull out DESeq2 results for multiple contrasts from httr_deg collection and reshape into wide matrix

Usage

```
getFCmatrix(
  db_host,
  db_name,
  collections = character(0),
  anl_name = getAnlName(mean_cnt = mean_cnt, plate_effect = plate_effect, shrink
      shrinkage),
  mean_cnt = getOption("mean_cnt"),
  plate_effect = T,
  shrinkage = "normal",
  stype = "test sample",
  debug = getOption("debug", default = FALSE),
  threads = -1,
  ...
)
```

Arguments

db_name	(character) = Specifies the DB name to connect to
collections	(named character vector) = Optional vector to re-map collection names used, where name=default collection; value=new collection name
anl_name	(<i>character</i>) = Which analysis configuration to pull results for, default: based on mean_cnt, plate_effect, and shrinkage params, but setting here takes precedence
mean_cnt	(integer) = passed to getAnlName to set anl_name, default: 5
plate_effect	(logical) = passed to getAnlName to set anl_name, default: T
shrinkage	(character) = passed to getAnlName to set anl_name, default: "normal"
stype	(<i>character</i>) = Filter for sample type, set to NULL to remove this filter, default: "test sample"
debug	(logical) = Whether to print debug messages, default: FALSE, overridden by options(debug=)
threads	(<i>integer</i>) = number of threads to use to parallelize = All other params are passed to mongoQuery to provide additional filtering (if chem_id is not specific enough, this will generate an error)
host	(character) = Specifies the host to open a connection to

Value

(*list of data.frame*) = First member \$fc is matrix of contrast (rows) x probe (col), Second member \$cmp is matrix of contrast (rows) x meta-data columns, both in same row order. Missing probes will be NA.

getProbeManifest 35

getProbeManifest

getProbeManifest Function to pull out all probe info as a data.frame This excludes the transcript field in order to get a flat table, but includes everything else The table is also sorted on the index field to ensure the row order of the original manifest

Description

getProbeManifest Function to pull out all probe info as a data.frame This excludes the transcript field in order to get a flat table, but includes everything else The table is also sorted on the index field to ensure the row order of the original manifest

Usage

```
getProbeManifest(
  DB = NULL,
  db_host = NULL,
  db_name = NULL,
  collection = "httr_probe",
  fields = "{\"transcripts\":0}",
  ...
)
```

Arguments

DB	(<i>mongo object</i>) = If specified, use this open connection to httr_probe collection, ignore db_host, db_name, collection
db_host	(character) = If DB is NULL, this specifies the host to open a connection to
db_name	(character) = If DB is NULL, this specifies the DB name to connect to
collection	(<i>character</i>) = If DB is NULL, this specifies the collection to connect to, default: httr_probe

Value

(data.frame) = Data frame with complete manifest

```
getQCdefault Function to determine appropriate threshold for a QC metric when not specified: NOTE: When the relevant flag is not in httrQCflags, this will always return NULL
```

Description

getQCdefault Function to determine appropriate threshold for a QC metric when not specified: NOTE: When the relevant flag is not in httrQCflags, this will always return NULL

```
getQCdefault(qc_metric)
```

36 getWellCounts

Arguments

```
qc_metric (character) = Name of the QC metric to get default for
```

Value

(numeric) = Value of default threshold, or NULL if not thresholding on this metric

getTrts

getTrts

• Low-level function to extract all docs from httr_trt_grp_cmp that match a specific query Before returning, the ctrl_wells and trt_wells members are converted to data frames

Description

getTrts

• Low-level function to extract all docs from httr_trt_grp_cmp that match a specific query Before returning, the ctrl_wells and trt_wells members are converted to data frames

Usage

```
getTrts(DB, query, debug = getOption("debug", default = FALSE))
```

Arguments

DB	(mongo object) = Open DB connection to httr_trt_grp_cmp collection
query	(json) = Query for returning docs from httr_trt_grp_cmp collection
debug	(logical) = Whether to print debug messages, default: FALSE, overridden by options(debug=)

Value

(*list of lists*) = Names come from trt_grp_id field, each member is a complete httr_trt_grp_cmp doc (list)

 ${\it getWellCounts}~{\it Function~to~pull~out~counts~and~treatment~info~as} \\ {\it data.frames~for~a~set~of~sample~IDs}$

Description

getWellCounts Function to pull out counts and treatment info as data.frames for a set of sample IDs

getWellInfo 37

Usage

```
getWellCounts(
  DB = NULL,
  db_host = NULL,
  db_name = NULL,
  collection = "httr_well",
  debug = getOption("debug", default = FALSE),
  sample_id = NULL,
  ...
)
```

Arguments

DB	(<i>mongo object</i>) = If specified, use this open connection to httr_well collection, ignore db_host, db_name, collection
db_host	(character) = If DB is NULL, this specifies the host to open a connection to
db_name	(character) = If DB is NULL, this specifies the DB name to connect to
collection	(<i>character</i>) = If DB is NULL, this specifies the collection to connect to, default: httr_well
debug	(<i>logical</i>) = Whether to print debug messages, default: FALSE, overridden by options(debug=) = Additional parameters passed to mongoQuery and used to specify which wells to extract counts for
sample_id	(<i>character</i>) = Used to query counts for a specific subset of sample IDs, if NULL this will not be part of query

Value

(*list*) = Has two members, both data.frames: \$treatments has rows = samples, \$counts has rows = probes, columns = samples

getWellInfo getWellInfo Function to pull out treatment info ONLY from httr_well as a data.frame for any relevant query

Description

getWellInfo Function to pull out treatment info ONLY from httr_well as a data.frame for any relevant query

```
getWellInfo(
  DB = NULL,
  db_host = NULL,
  db_name = NULL,
  collection = "httr_well",
  debug = getOption("debug", default = FALSE),
  ...
)
```

38 get_global_option

Arguments

DB	(mongo object) = If specified, use this open connection to httr_well collection, ignore db_host, db_name, collection
db_host	(character) = If DB is NULL, this specifies the host to open a connection to
db_name	(character) = If DB is NULL, this specifies the DB name to connect to
collection	$(character) = If DB$ is NULL, this specifies the collection to connect to, default: httr_well
debug	(<i>logical</i>) = Whether to print debug messages, default: FALSE, overridden by options(debug=) = Any additional parameters are passed to mongoQuery to constrain the query

Value

(*data.frame*) = Table of well treatment info from httr_well, currently excludes _id, probe_cnts, raw_id, counts_id, trt_id, sorts by sample_id

 $\verb|get_global_option| \textit{ get_global_option Function to find the value of the passed global_var}|$

first checks if old name for global var exists and returns that and a warning if not, then just simply return the value corresponding to the

passed global var

Description

get_global_option Function to find the value of the passed global_var first checks if old name for global var exists and returns that and a warning if not, then just simply return the value corresponding to the passed global var

Usage

```
get_global_option(global_var)
```

Arguments

global_var (character) = may be any of the following values: min_mapd_frac, min_n_reads_mapd, min_n_sig80, min_n_cov5, max_top10_prop, max_gini_coef, qc_flags

Value

numeric

insertByID 39

insertByID	title insertByID - Insert a list of documents into a MongoDB collection without creating redundant entries. Takes a list of dicts and first checks whether any documents with matching IDs are already present.
	If so, these documents are either skipped for insert (rerun=False) or replaced (rerun=True)

Description

title insertByID - Insert a list of documents into a MongoDB collection without creating redundant entries. Takes a list of dicts and first checks whether any documents with matching IDs are already present. If so, these documents are either skipped for insert (rerun=False) or replaced (rerun=True)

Usage

```
insertByID(
   DB,
   docs,
   id_field = "sample_id",
   rerun = FALSE,
   debug = getOption("debug", default = FALSE)
)
```

Arguments

```
DB (mongo object) = Open database/collection connection for inserting

docs (list) = List of documents to insert into collection

id_field (character) = The field to use for matching up documents, default: sample_id

rerun (logical) = Whether to overwrite existing data with same sample_id

debug (logical) = Whether to report debug messages
```

Value

(integer) = The number of documents successfully inserted, Note: Unlike the pymongo API, mongolite does not return the IDs after a successful insert NOTE: For this function in particular, it makes sense to have a custom wrapper for each collection, e.g. to make sure there are matching IDs in other collections

```
insertManyDEG insert multiple docs into httr_deg collection, matched by combination of trt_grp_id and anl_name will skip docs that are already inserted if rerun = False
```

Description

insertManyDEG Insert multiple docs into httr_deg collection, matched by combination of trt_grp_id and anl_name will skip docs that are already inserted if rerun = False

40 insertOneDEG

Usage

```
insertManyDEG(
  DB = NULL,
  deg_docs = NULL,
  db_host = NULL,
  db_name = NULL,
  collection = "httr_deg",
  rerun = FALSE,
  debug = getOption("debug", default = FALSE),
  ...
)
```

Arguments

DB	(<i>mongo object</i>) = If specified, use this open connection to httr_deg collection, ignore db_host, db_name, collection
deg_docs	$(list\ of\ lists)$ = Each top-level entry should be a document (list) corresponding to httr_deg schema
db_host	(emphcharacter) = If DB is $NULL$, this specifies the host to open a connection to
db_name	(character) = If DB is NULL, this specifies the DB name to connect to
collection	(character) = If DB is NULL, this specifies the collection to connect to, default: httr_deg
rerun	(<i>logical</i>) = Whether to re-run the insertion of this data and override existing doc with matching trt_grp_id, anl_name
debug	(logical) = Whether to print debug messages, default: FALSE, overridden by options(debug=) = Additional params passed to insertOneDEG

Value

(integer) = Number of documents successfully inserted

```
insertOneDEG insertOneDEG - Insert a doc into httr_deg collection, matched by combination of trt_grp_id and anl_name
```

Description

insertOneDEG - Insert a doc into httr_deg collection, matched by combination of trt_grp_id and anl_name

```
insertOneDEG(
  DB = NULL,
  deg_doc = NULL,
  db_host = NULL,
  db_name = NULL,
  collection = "httr_deg",
```

insert_into_httr_study 41

```
init_note = "initial run",
  update_note = NULL,
  rerun = FALSE,
  debug = getOption("debug", default = FALSE)
```

Arguments

DB	(mongo object) = If specified, use this open connection to httr_deg collection, ignore db_host, db_name, collection
deg_doc	(list) = Single document corresponding to httr_deg schema for insert
db_host	(character) = If DB is NULL, this specifies the host to open a connection to
db_name	(character) = If DB is NULL, this specifies the DB name to connect to
collection	$(character)$ = If DB is NULL, this specifies the collection to connect to, default: httr_deg
init_note	(character) = If this is first insert into DB, set update_notes.note to this value
update_note	$(character)$ = If this is subsequent insert into DB, set last update_notes.note to this value (defaults to: "update on MM/DD/YY")
rerun	(<i>logical</i>) = Whether to re-run the insertion of this data and override existing doc with matching trt_grp_id, anl_name
debug	(logical) = Whether to print debug messages, default: FALSE, overridden by options(debug=)

Value

(logical) = TRUE if inserted successfully, FALSE otherwise

Description

insert_into_httr_study This function reads all inputs, validates them to create a document into httr_study

```
insert_into_httr_study(
   study_id = "httr_mcf7_pilot",
   study_name = "HTTr MCF-7 Pilot Study",
   study_desc =
     "Pilot of HTTr platform on MCF-7 cells testing multiple media types and expostudy_probe = "httr_probe",
   study_well = "httr_well",
   study_well = "httr_well",
   study_trt_grp_cmp = "httr_trt_grp_cmp",
   study_degs = "httr_deg",
   src = "BioSpyder",
```

```
tech = "TempoSeq",
 times_hr = list(6, 12, 24),
 max\_dose\_level = 8,
 media = list("DMEM", "PRF.DMEM"),
 cell_types = list("MCF-7"),
 chem_ids = list("TP0001651A01", "TP0001651A02", "TP0001651A03", "TP0001651A04"
    "TP0001651A05", "TP0001651A06", "TP0001651B01", "TP0001651B02", "TP0001651B0
    "TP0001651B04", "TP0001651B05", "TP0001651B06", "TP0001651C01", "TP0001651C0
    "TP0001651C03", "TP0001651C04", "TP0001651C05", "TP0001651C06"),
 nreps = 3,
 bs_assay_id = "",
 assay_name = "human_wt_1.2_mcf7",
 platform name = "human wt",
 platform\_ver = "1.2",
 atten_type = "mcf7",
 probe_source_file = "180905 Human Whole Transcriptome 1.2 Manifest.xlsx",
 probe_source_date = "2018-09-05T00:00:00Z",
 probe_source_rcv = "2019-07-15T13:24:12Z",
 probe_ncct_file = "httr_mcf7_pilot_probe.csv",
 probe_ncct_date = "2019-09-23T13:44:50.419Z",
  fasta file =
    "/share/projects/HTTr/HTTr_pipeline_dev/HumanWT_v1/AlignCount/genome/humanWT
 anno_type = "refseq",
 anno_date = "2018-09-05T00:00:00Z",
 host = NULL
 db = NULL
)
```

Arguments

```
(str) = short string with study ID (all collections specific to this study should
study_id
                  have name starting with this ID, so no spaces, punctuation or special characters)
                  (str) = Name of the study in human readable format (spaces, punctuation, etc.
study_name
                  allowed) – replaces Imran's "name" field in prev version of this collection.
                  (str) = longer text briefly describing the study
study_desc
study_probe
                  (str) = name of collection that catalogs all probes in the version of probe set used
                  for this study
                  (str) = name of collection with the well-level data for this study
study_well
study_trt_grp_cmp
                  (str) = name of collection with the treatment groups for comparison for this
                  (str) = name of collection with the DEG data for this study
study degs
                  (str) = Source of data, e.g. "BioSpyder"
src
                  (str) = Technology type, e.g. "TempoSeq"
tech
                  (list of int) = List of exposure times used (does not include Bulk Lysates or other
times hr
                  QC samples)
max_dose_level
```

(float) = Highest dose level for any conc-response curve - this replaces Imran's concs_um field because not all concs are um or will match up across chemicals, but knowing the intended number of points in each conc-response curve is useful.

insert_into_httr_study 43

media	(list of str) = List of media tested (does not include Bulk Lysates or other QC samples
cell_types	(list of str) = List of cell-types tested (was just "cell" in previous schema; does not include Bulk Lysates or other QC samples)
chem_ids	(list of str) = List of chemical IDs tested (whatever IDs are used in chem_id field of httr_well_trt collection).
nreps	(int) = Number of replicates for primary treatment groups of interest (for screening studies, this is per concentration of test samples, also per media/time/cell-type in other study types)
bs_assay_id	(str) = Unique bar code or other identifier for the actual tube/batch of assay mix that BioSpyder used - we will need to get this information from them for tracking
assay_name	(str) = Unique identifier for this probe set version (including the combination of attenuation factors, e.g. human_wt_1.2_mcf7 to denote Human Whole Transcriptome v1.2 with MCF-7 attenuation)
platform_nam	
	(str) = Just the general platform, without version number or attenuation modifiers, e.g. human_wt for human whole transcriptome, or human_s1500
platform_ver	(str) = Version number for the platform, e.g. 1.2 or 2.0
atten_type	(str) = Type of attenuation used, e.g. "none" if no attenuation, "mcf7", "heparg_super", or "broad"
probe_source	
	(str) = Name of the manifest file sent from BioSpyder that was used to process this data
probe_source	
1	(date) = When did BioSpyder last update the file?
probe_source	
nroho naat f	(date) = When did we actually receive the manifest from BioSpyder?
probe_ncct_f	(str) = If NCCT performed additional clean-up of the manifest, name of the
	final clean file that was loaded into DB. Otherwise can be blank or same as source_file.
probe_ncct_d	
	(date) = When did NCCT finalize the probe set manifest (after validating/modifying the manifest received from BioSpyder)
fasta_file	(str) = Full path to fasta file used to build HISAT2 index and align fastq files
anno_type	(str) = What was the primary annotation used for this probe set (e.g. earlier versions are refseq, newer versions are ensembl)
anno_date	(date) = When was the annotation for this probe set last updated (we may update the annotations periodically to keep up with current transcriptome annotations)
	process: The function does validation checks on an entry into the httr_study

Value

nothing explicit print the test that fail if any. If validation passes, nothing returned

44 mongoQuery

iterate	iterate with DB object as well as fields and potentially query named
	list, inquire db for matching documents

Description

iterate with DB object as well as fields and potentially query named list, inquire db for matching documents

Usage

```
iterate(DB = NULL, query = NULL, fields = NULL)
```

Arguments

DB (mongo object) = Opened connection to specific Mongo Database and Collec-

tion

fields: (json list) of fields to be included in our search

query: (json query) to be used in our search

Value

list of documents matching search specified by query/fields

mongoQuery Build a mongo query from an arbitrary set of params or a list Takes a list of named values/vectors and converts to a JSON-

a list Takes a list of named values/vectors and converts to a JSONformat mongo query. Each entry with length > 1 is interpreted as a set of potential matching values and converted to the \$in keyword for

mongo

Description

mongoQuery Build a mongo query from an arbitrary set of params or a list Takes a list of named values/vectors and converts to a JSON-format mongo query. Each entry with length > 1 is interpreted as a set of potential matching values and converted to the \$in keyword for mongo

Usage

```
mongoQuery(...)
```

Arguments

... (any) = Any set of named parameters OR a single named list of all arguments

mongoURL 45

mongoURL MongoURL

Description

Just format a mongo URL based on mongoServer, username, password, and database

Usage

```
mongoURL(host, user, passwd, db, authSource = NULL, authMechanism = NULL)
```

Arguments

```
host (character)
user (character)
passwd (character)
db (character)
authSource (character)
authMechanism (character)
```

Value

formatted string to pass to mongolite library connect function

openMongo

openMongo open a connection to a specific collection in a mongo db

Description

openMongo open a connection to a specific collection in a mongo db

```
openMongo(
  host = getOption("httrDefaultHost"),
  user = NULL,
  passwd = NULL,
  db = NULL,
  collection = NULL,
  authSource = NULL,
  authMechanism = NULL
)
```

46 outerFences

Arguments

host	(character)
user	(character)
passwd	(character)
db	(character)
collection	(character)
authSource	(character)
authMechanis	sm
	(character)

Value

formatted string to pass to mongolite library connect function

outerFences	outerFences Function to compute "Tukey's Outer Fence" for a set of
	values Computes both lower and upper thresholds, which are 3*IQR
	beyond the lower and upper quartile respectively

Description

outerFences Function to compute "Tukey's Outer Fence" for a set of values Computes both lower and upper thresholds, which are 3*IQR beyond the lower and upper quartile respectively

Usage

```
outerFences(x, iqr.factor = 3)
```

Arguments

```
x (numeric vector) = Set of values to compute fences for iqr.factor (numeric) = Multiplier for IQR in the fence computation, default is 3
```

```
(numeric, length 2) = Lower and upper fence values
```

plotQCjitter 47

plotQCjitter	plotQCjitter Function for drawing jitter plot of a QC metric grouped by an experimental design factor

Description

plotQCjitter Function for drawing jitter plot of a QC metric grouped by an experimental design factor

Usage

```
plotQCjitter(
   qc_metric,
   design_fac,
   data,
   cutoff = getQCdefault(qc_metric),
   alt_cutoff = NULL,
   title = NULL,
   plot.ylim = NULL
)
```

Arguments

qc_metric	(character) = Name of the QC metric column in data to plot
design_fac	(character) = Name of the design factor column in data to group by
data	(data.frame) = Table of all QC/meta-data for plotting
cutoff	(<i>numeric</i>) = Single value for drawing a horizontal dashed line, denoting current threshold - if NULL, will attempt to fill in with appropriate default
alt_cutoff	(numeric) = Single value for drawing a horizontal dotted line, denoting alternate threshold - if NULL, not drawn
title	(character) = Title for the plot, optional
plot.ylim	(<i>numeric</i> , length 2) = Set lower and upper limits for Y axis for consistent plotting, if NULL plots will auto-decide

Value

(ggplot2) object that can be further modified or drawn with print

plotQCviolin Function for drawing violin plot of a QC metric grouped by an experimental design factor
oy an experimental aesign factor

Description

plotQCviolin Function for drawing violin plot of a QC metric grouped by an experimental design factor

48 process_chemInfo

Usage

```
plotQCviolin(
   qc_metric,
   design_fac,
   data,
   cutoff = getQCdefault(qc_metric),
   alt_cutoff = NULL,
   title = NULL,
   plot.ylim = NULL
)
```

Arguments

qc_metric	(character) = Name of the QC metric column in data to plot
design_fac	(character) = Name of the design factor column in data to group by
data	(data.frame) = Table of all QC/meta-data for plotting
cutoff	(numeric) = Single value for drawing a horizontal dashed line, denoting current threshold - if NULL, will attempt to fill in with appropriate default
alt_cutoff	(numeric) = Single value for drawing a horizontal dotted line, denoting alternate threshold - if NULL, not drawn
title	(character) = Title for the plot, optional
plot.ylim	(<i>numeric</i> , length 2) = Set lower and upper limits for Y axis for consistent plotting, if NULL plots will auto-decide

Value

(ggplot2) object that can be further modified or drawn with print

```
process_chemInfo process_chemInfo
```

Description

```
process_chemInfo
```

Usage

```
process_chemInfo(wellTrt)
```

Arguments

```
wellTrt: (character)data.table with data to check
```

Value

chemInfo data table a subset of the inputted data.table to be stored into httr_chem by default

```
process_chem_columns
```

process_chem_columns takes as input: (coming from prior manipulations)

Description

process_chem_columns takes as input: (coming from prior manipulations)

Usage

```
process_chem_columns(wellTrt)
```

Arguments

wellTrt:

(character)data.table with data to check

Value

modified data table after removing rows with dup chems and removing chem_name, dtxsid, casrn and bottle_id columns

qcBatch

qcBatch Function to run QC functions on all documents in httr_counts and store results in httr_counts_qc Also pulls probe info out of httr_probe (should delegate to functions in db/probes.R), NOTE: There's no real reason here to do sample_id batching, but could still use config files to store the DB info and collection names?

Description

qcBatch Function to run QC functions on all documents in httr_counts and store results in httr_counts_qc Also pulls probe info out of httr_probe (should delegate to functions in db/probes.R), NOTE: There's no real reason here to do sample_id batching, but could still use config files to store the DB info and collection names?

```
qcBatch(
  db_host,
  db_name,
  db_probe = "httr_probe",
  db_counts = "httr_counts",
  db_collection = "httr_counts_qc",
  keep_flags = "OK",
  calc_flag = TRUE,
  rerun = FALSE,
  min_mapd_frac = get_global_option("min_mapd_frac"),
  min_n_reads_mapd = get_global_option("min_n_reads_mapd"),
```

50 readWellTrtFile

```
min_n_sig80 = get_global_option("min_n_sig80"),
min_n_cov5 = get_global_option("min_n_cov5"),
max_top10_prop = get_global_option("max_top10_prop"),
max_gini_coef = get_global_option("max_gini_coef"),
qc_flags = get_global_option("qc_flags"),
debug = getOption("debug", default = FALSE),
...
)
```

Arguments

```
(character) = MongoDB host URL
db_host
db_name
                 (character) = MongoDB name
                 (character) = Name of probe collection, default: httr_probe
db_probe
                 (character) = Name of counts collection, default: httr counts
db counts
db collection
                 (character) = Name of counts_qc collection, default: httr_counts_qc
                 (character vector) = Keep all probes with these flags in the downstream QC
keep_flags
                 (everything else is considered BAD), default: "OK"
                 (logical) = Whether or not to calculate qc_flag, default True but when set to
calc_flag
                 False all flags will be set to "OK"
                 (logical) = Whether to overwrite existing data in httr_counts_qc
rerun
min_mapd_frac
                 (numeric) = minimum mapd_frac, flag with LOW_MAPD_FRAC below this
min_n_reads_mapd
                 (integer) = minimum n_reads_mapd, flag with LOW_READS below this cutoff
                 (integer) = minimum n_sig80, flag with LOW_NSIG80 below this cutoff
min_n_sig80
min n cov5
                 (integer) = minimum n cov5, flag with LOW NCOV5 below this cutoff
max_top10_prop
                 (numeric) = maximum top10_prop, flag with HIGH_TOP10 above this cutoff
max_gini_coef
                 (numeric) = maximum gini_coef, flag with HIGH_GINI above this cutoff in the
                 ... param, a query can be passed that will select the documents in httr_well_trt
                 used to further filter the sample_ids to run the qc_flag calculation and storage in
                 the httr_counts_qc collection
                 (logical) = Whether to report debug messages
debug
```

readWellTrtFile readWellTrtFile takes as input: (coming from prior manipulations)

Description

readWellTrtFile takes as input: (coming from prior manipulations)

```
readWellTrtFile(sampleID_file)
```

refChemGroup 51

Arguments

Value

data table with json data with corrected column type

refChemGroup Construct a document for httr_trt_grp_cmp collection corresponding to a specific reference chemical treatment.

Description

Given a specific reference chemical and plate group (or media/timeh that specify plate group), plus any optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This function is just a convenenience wrapper to chemTreatGroup with modified defaults.

Usage

```
refChemGroup(
   DB,
   trt_chem,
   pg_id = NULL,
   media = NULL,
   timeh = NULL,
   trt_dose = NULL,
   trt_type = "reference chemical",
   trt_prop_fields = c("chem_id", "conc", "conc_unit", "stype"),
   ...
)
```

Arguments

```
DB
                   (string) = name of database with httr_well and httr_trt_grp_cmp collections col-
                   lections
                   (str) = Match to chem_id in httr_well
trt_chem
pg_id
                   (str) = Filter both trt and ctrl wells, match to pg_id in httr_well, defaults to None
media
                   (str) = Filter both trt and ctrl wells, match to media in httr_well, defaults to None
                   (str) = Filter both trt and ctrl wells, match to timeh in httr_well, defaults to None
timeh
trt_dose
                   (int) = Match to dose_level in httr_well - Excluded by default
                   (str) = Match to stype in httr_well, defaults to "reference chemical"
trt_type
trt_prop_fields
```

(vector) = Fields in httr_well collection that should match up for all trt_wells and the singular value should be propagated to field of same name in httr_trt_grp_cmp collection. Default is to propagate chem_id, conc, conc_unit, and stype fields (exclude dose_level here). ... = All additional args passed to treatGroupFrom-Samples

52 refRNAGroup

Details

Parameters:

Value

(list) = Document that was inserted into trt_grp collection Warn if neither pg_id nor media+timeh were specified:

refRNAGroup

refRNAGroup Construct a document for httr_trt_grp_cmp collection corresponding to a specific reference RNA QC sample comparison.

Description

Given a specific plate group, plus any optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This function is just a convenenience wrapper to chemTreatGroup with modified defaults.

Usage

```
refRNAGroup(
   DB,
   pg_id,
   trt_type = "QC sample",
   ctrl_type = "QC sample",
   trt_src = "HBRR",
   ctrl_src = "UHRR",
   trt_prop_fields = c("stype"),
   both_prop_fields = c("pg_id", "block_id"),
   ctrl_desc_field = "trt_name",
   grp_id_opts = c("vs", "pg"),
   ...
)
```

Arguments

```
DB (string) = name of database with httr_well and httr_trt_grp_cmp collections collections collections collections collections (str) = Filter both trt and ctrl wells, match to pg_id in httr_well, required for ref RNA groups

trt_type (str) = Match to stype in httr_well, defaults to "QC sample"

ctrl_type (str) = Match to stype in httr_well, defaults to "QC sample"

trt_src (str) = Match to rna_src in httr_well, defaults to "HBRR"

ctrl_src (str) = Match to rna_src in httr_well, defaults to "UHRR"

trt_prop_fields
```

(list) = Fields in httr_well collection that should match up for all trt_wells and the singular value should be propagated to field of same name in httr_trt_grp_cmp collection. Default is to propagate just the stype field (exclude all chemical and dose information).

runDESeq2 53

both_prop_fields

(list of str) = Fields in httr_well collection that should match up for all trt_wells AND ctrl_wells, and the singular value should be propagated to to field of same name in httr_trt_grp_cmp collection. Default is to propagate pg_id and block_id only for ref RNA comparisons

ctrl_desc_field

(str) = Field to propagate from ctrl_wells to ctrl field in httr_trt_grp_cmp document, defaults to "trt_name" for ref RNA comparisons

grp_id_opts

(list) = Options to modify trt_grp_id generation, "vs" will capture the pairwise comparison type, and "pg" is usually sufficient for ref RNA to create distinct group IDs for each plate group ... = All additional args passed to treat-GroupFromSamples

Details

Parameters:

Value

(list) = Document that was inserted into trt_grp collection """ Pass to chemTreatGroup:

runDESeq2

runDESeq2 - Core function for running DESeq2 on a set of wells This was adapted from Imran's previous DESeq2 function and I tried to keep all parameters the same Previous version can be found in httr-wf-dev repo in BitBucket, httr-ph-i branch, lib/gexp/deseq2.py, run-DESeq4 function The main change here is that the default shrinkage method is now done with the lfcShrink function in DESeq2 package, instead of the old method The column used for primary treatment groups is now trt_name instead of trt_id to keep with the DB schema NOTE: This function does no sample or probe filtering - that was implemented on the python side in Imran's code, and will be handled by a separate function in the R lib NOTE: Derik's version converted conc as a factor and used that as the primary model - results should be equivalent but might be a good idea to sort samples in increasing conc order, ideally in a higher level function that's specific to handling dose-response sets

Description

runDESeq2 - Core function for running DESeq2 on a set of wells This was adapted from Imran's previous DESeq2 function and I tried to keep all parameters the same Previous version can be found in httr-wf-dev repo in BitBucket, httr-ph-i branch, lib/gexp/deseq2.py, runDESeq4 function The main change here is that the default shrinkage method is now done with the lfcShrink function in DESeq2 package, instead of the old method The column used for primary treatment groups is now trt_name instead of trt_id to keep with the DB schema NOTE: This function does no sample or probe filtering - that was implemented on the python side in Imran's code, and will be handled by a separate function in the R lib NOTE: Derik's version converted cone as a factor and used that as the primary model - results should be equivalent but might be a good idea to sort samples in increasing cone order, ideally in a higher level function that's specific to handling dose-response sets

Usage

```
runDESeq2(
   COUNTS,
   CONDS,
   ref_level = "DMSO_0",
   plate_effect = F,
   plate_cont = F,
   shrinkage = "normal",
   threads = 1,
   debug = getOption("debug", default = FALSE)
)
```

Arguments

COUNTS	(<i>data.frame</i>) = Data frame containing all counts to be analyzed, row.names = probes, colnames = samples
CONDS	(<i>data.frame</i>) = Data frame containing treatment data for all samples, rows = samples in same order as COUNTS, must have trt_name and plate_id columns
ref_level	(character) = Value of "reference" treatment in trt_name column (will be used as model intercept)
plate_effect	(logical) = Should the DESeq2 model include plate effects?
plate_cont	(logical) = Should the plate contrasts be returned in addition to the primary treatment contrasts? (CURRENTLY NOT SUPPORTED)
shrinkage	(character) = Type of shrinkage to use, current options are "normal" for new default shrinkage, "none" for no shrinkage
threads	(integer) = Number of threads to use, if > 1 parallelization is accomplished with BiocParallel package
debug	(logical) = Whether to print debug messages, default: FALSE, overridden by options(debug=)

Value

(*list of data.frame*) = List with two members, first is the data.frame of all results for primary treatments, second is for plate effects Each data.frame has columns probe_id, trt_name, baseMean, log2FoldChange, lfcSE, stat, pvalue, padj

 $\verb"runDESeq2ForChemCond"$

runDESeq2ForChemCond - Run DESeq2 on a single chemical doseresponse series Uses DB access functions to extract all httr_trt_grp_id entries corresponding to a dose-response series, then extracts counts for all wells, filters by min read depth and mean floor, and passes through runDESeq2. Also optionally saves the results back to the DB in addition to returning a single data.frame with all results that can be written to a TSV file or saved in Rdata format.

Description

runDESeq2ForChemCond - Run DESeq2 on a single chemical dose-response series Uses DB access functions to extract all httr_trt_grp_id entries corresponding to a dose-response series, then extracts counts for all wells, filters by min read depth and mean floor, and passes through runDESeq2. Also optionally saves the results back to the DB in addition to returning a single data.frame with all results that can be written to a TSV file or saved in Rdata format.

Usage

```
runDESeq2ForChemCond(
  chem_id,
  db_host,
  db_name,
 pg_id = NULL,
 media = NULL,
  timeh = NULL,
 max_dose_level = 8,
 min colsum = 10^5,
  mean cnt = getOption("mean cnt"),
 plate_effect = F,
  shrinkage = "normal",
  threads = 1,
  collections = character(0),
  db_insert = T,
  rerun = F,
 note = NULL,
  debug = getOption("debug", default = FALSE)
)
```

Arguments

```
(character) = ID for the chemical to run, matched against chem_id field in
chem_id
                  httr_trt_grp_cmp
db_host
                  (character) = Host for DB connection
                  (character) = Name of DB
db_name
pg_id
                  (character) = Limit httr_trt_grp_cmp to this plate group, useful if multiple dose
                  response series for same chem id
media
                  (character) = Limit httr_trt_grp_cmp to this media type, useful if multiple dose
                  response series for same chem_id
                  (integer) = Limit httr_trt_grp_cmp to this timeh, useful if multiple dose response
timeh
                  series for same chem_id
max_dose_level
                  (integer) = Max number of dose levels expected
                  (integer) = Filter samples by min read total, set to NULL to skip, defaults to
min_colsum
                  100k
mean_cnt
                  (integer) = Filter probes by mean floor, set to NULL to skip, defaults to 5
plate_effect (logical) = Whether or not to include plate effect in model, defaults to FALSE
                  (character) = Type of shrinkage to use, defaults to new default, a.k.a. "normal"
shrinkage
                  (integer) = Number of threads to use for core DESeq2 functions, defaults to 1
threads
                  (no multi-threading)
```

56 runDESeq2Single

collections	(named character vector) = Optional vector to re-map collection names used, where name=default collection; value=new collection name
db_insert	(logical) = Whether to insert the results back into the DB, defaults to TRUE
rerun	(logical) = Whether to re-run the DB insertion if results are in there already (ideally this should check the DB before running)
note	(character) = Overwrite default note fields for insert if this is specified.
debug	(<i>logical</i>) = Whether to print debug messages, default: FALSE, overridden by options(debug=)

Value

(data.frame) = Contains DESeq2 results for all treatments in the dose-response series

runDESeq2Single

runDESeq2Single - Run DESeq2 on a single contrast (e.g. single conc reference chemical, bulk lysate, or ref RNA QC samples) Uses DB access functions to extract the httr_trt_grp_cmp entry corresponding to a specific comparison, then extracts counts for all wells, filters by min read depth and mean floor, and passes through runDESeq2. Also optionally saves the results back to the DB in addition to returning a single data.frame with all results that can be written to a TSV file or saved in Rdata format.

Description

runDESeq2Single - Run DESeq2 on a single contrast (e.g. single conc reference chemical, bulk lysate, or ref RNA QC samples) Uses DB access functions to extract the httr_trt_grp_cmp entry corresponding to a specific comparison, then extracts counts for all wells, filters by min read depth and mean floor, and passes through runDESeq2. Also optionally saves the results back to the DB in addition to returning a single data frame with all results that can be written to a TSV file or saved in Rdata format.

```
runDESeq2Single(
  trt_grp_id,
  db_host,
  db_name,
  min_colsum = 10^5,
  mean_cnt = getOption("mean_cnt"),
  plate_effect = F,
  shrinkage = "normal",
  threads = 1,
  collections = character(0),
  db_insert = T,
  rerun = F,
  note = NULL,
  debug = getOption("debug", default = FALSE),
  ...
)
```

sameCtrlWells 57

Arguments

trt_grp_id	(character) = ID for the trt_grp_cmp to run, matched against trt_grp_id field in httr_trt_grp_cmp
db_host	(character) = Host for DB connection
db_name	(character) = Name of DB
min_colsum	(integer) = Filter samples by min read total, set to NULL to skip, defaults to $100k$
mean_cnt	(integer) = Filter probes by mean floor, set to NULL to skip, defaults to 5
plate_effect	(logical) = Whether or not to include plate effect in model, defaults to FALSE
shrinkage	(character) = Type of shrinkage to use, defaults to new default, a.k.a. "normal"
threads	(integer) = Number of threads to use for core DESeq2 functions, defaults to 1 (no multi-threading)
collections	(named character vector) = Optional vector to re-map collection names used, where name=default collection; value=new collection name
db_insert	(logical) = Whether to insert the results back into the DB, defaults to TRUE
rerun	(logical) = Whether to re-run the DB insertion if results are in there already (ideally this should check the DB before running)
note	(character) = Overwrite default notes with this one in update_notes field.
debug	$(logical) = \mbox{Whether to print debug messages, default: FALSE, overridden by options(debug=)}$

Value

(data.frame) = Contains DESeq2 results for all treatments in the dose-response series

sameCtrlWells sameCtrlWells

• Function to check that a list of httr_trt_grp_cmp docs all have the same control group This is useful when extracting treatment groups that should correspond to a dose-response series

Description

sameCtrlWells

• Function to check that a list of httr_trt_grp_cmp docs all have the same control group This is useful when extracting treatment groups that should correspond to a dose-response series

Usage

```
sameCtrlWells(trt_list, debug = getOption("debug", default = FALSE))
```

Arguments

Value

(logical) = FALSE if there are multiple entries in trt_list with different ctrl_well sets - automatically TRUE if length(trt_list) < 2 although this will trigger a debug message

```
sampleID_scan_and_update

sampleID_scan_and_update This function takes two files and scan
them for differences
```

Description

sampleID_scan_and_update This function takes two files and scan them for differences

Usage

```
sampleID_scan_and_update(
  orig_sampleID_file = NULL,
  sampleID_file = NULL,
  targetWellTrtCol = "httr_well_trt",
  targetChem = "httr_chem",
  skipped_tests = skipped_tests,
  db_host = host,
  db_name = db,
  ...
)
```

Arguments

```
orig_sampleID_file:
                  (character) file that has already been processed prior
sampleID_file:
                  (character) new vs of that files with some edits - regardless where the edits took
                  place
targetWellTrtCol:
                  (character) target collection where the updates will be written to - by default
                  httr_well_trt
                  (character) the target Chem collection where some updates will be written to -
targetChem:
                  by default httr_chem
skipped_tests:
                  (character vector) the list of tests to skip when validating the resulting set of
db_host:
                  (character) mongo db server +- port
db_name:
                  (character) mongo db database
                  process: The function does a diff on the two files and calls sampleID_wrapper()
                  on that resulting list of modified sampleIds if the list validates against a list of
                  provided tests, the function continues to replace found sampleId documents in
                  target collections.
```

Value

nothing explicit print the test that fail if any. If validation passes, nothing returned

sampleID_wrapper 59

```
{\tt sampleID\_wrapper} \quad \textit{sampleID\_wrapper takes either as input: (coming from prior manipulations)}
```

Description

sampleID_wrapper takes either as input: (coming from prior manipulations)

Usage

```
sampleID_wrapper(
  wellTrt = NULL,
  chemInfo = NULL,
  sampleID_file = NULL,
  targetWellTrtCol = "httr_well_trt",
  targetChem = "httr_chem",
  status = "KEEP",
  skipped_tests = c(),
  db_host = NULL,
  db_name = NULL,
  validate = TRUE,
  max_dose_level = 8,
  required_cols = c("sample_id", "plate_id", "well_id", "trt_name", "qc_flag"),
  extra_cols = c(),
  ...
)
```

Arguments

Target

Chem Collection: (character) targetChem defaulting to httr_chem

status: (numeric) whether of not we keep prior data if found in target collection (sta-

tus="KEEP"), erase all items in collection (status="RERUN"), or just replace

the sampleIDs found (status="REPLACE")

tests: (character vector) list of tests that will be performed in the validation step, an

example is provided in tests.json

db_host: (character) mongo url with or without port

db_name: (character) mongo database or sandbox

60 tests

Value

nothing explicit

the function will process each sample ids, first making sure they pass the validation tests if they don't pass any of the validation tests declared in the relevant_test list, the offended test will be printed out and the the process is terminated if the validation tests pass, the process continues on to store (by erasing first all or some or none of potentially already present sampleIDs in the target collections) the data into two collections: the target collection and the target Chem collection as specified earlier finally, each sampleID that ended up being replaced in the target collection is removed from the httr_well collection, so that another downstream process may detect the dicrepencies between httr_well_trt and httr_well and rebuild the missing sampleIDs that were found in httr_well_trt but missing in httr_well

splitIDs Title splitIDs - Split a list of IDs based on which are present/absent in a collection. Given a list of IDs, search against a particular collection, and return a dict with two lists keyed by "present" and "absent".

Description

Title splitIDs - Split a list of IDs based on which are present/absent in a collection. Given a list of IDs, search against a particular collection, and return a dict with two lists keyed by "present" and "absent".

Usage

```
splitIDs(DB, ids, id_field = "sample_id")
```

Arguments

DB (mongo object) = Open connection to a Mongo DB and Collection

ids (vector) = Vector of IDs (typically character)

id_field (character) = Name of ID field, default: sample_id

Value

(*list*) with two members: \$present = vector of IDs that were found in the collection \$absent = vector of IDs that were NOT found in collection

tests validate httr well trt schema

Description

```
validate_httr_well_trt_schema
```

Usage

tests

Arguments

wellTrt:

(character vector) a list of named values, such as:

sample_id plate_id well_id chem_id chem_name dtxsid casrn conc conc_unit dose_level block_id replicate_num pg_id culture_id cell_type media stype trt_name rna_src qc_flag TC00001203_C07 TC00001203 C07 EPAPLT0593A13222 Glybenclamide DTXSID0037237 10238-21-8 100.00 uM 8 1 1 1 c2021-03-15 U-2 OS DMEM + 10% FBS test sample EPAPLT0593A13_8_100uM Plated Cells OK

skipped_tests:

(character vector) a list of skipped tests as defined in the tests.json file that are not used to test validate against

this long function performs many tests on the data quality and integrity provided in the wellTrt list it delineates the tests that have failed

Format

An object of class list of length 51.

Value

nothing explicit

treatGroupFromSamples

treatGroupFromSamples Construct a document for httr_trt_grp_cmp collection. Given a list of treatment and control samples, construct an appropriate document and insert into httr_trt_grp_cmp Parameters:

Description

treatGroupFromSamples Construct a document for httr_trt_grp_cmp collection. Given a list of treatment and control samples, construct an appropriate document and insert into httr_trt_grp_cmp Parameters:

```
treatGroupFromSamples(
   DB = getOption("DB_NAME"),
   trt_wells,
   ctrl_wells,
   trt_prop_fields = c("chem_id", "conc", "conc_unit", "dose_level", "stype"),
   both_prop_fields = c("media", "timeh", "pg_id", "block_id"),
   ctrl_desc_field = "stype",
   grp_id_opts = c(),
   well = "httr_well",
   trt_grp = "httr_trt_grp_cmp",
   rerun = FALSE,
   db_insert = TRUE,
   ...
)
```

62 treatPairGroup

Arguments

(string) = name of database with httr_well and httr_trt_grp_cmp collections DB (vector) = Sample IDs of all treated wells in the comparison group trt_wells (vector) = Sample IDs of all control wells in the the comparison group ctrl wells trt_prop_fields (vector) = Fields in httr_well collection that should match up for all trt_wells and the singular value should be propagated to field of same name in httr_trt_grp_cmp collection. Default is to propagate chem_id, conc, conc_unit, dose_level, and stype fields. both_prop_fields (vector) = Fields in httr_well collection that should match up for all trt_wells AND ctrl_wells, and the singular value should be propagated to to field of same name in httr_trt_grp_cmp collection. Default is to propagate media, timeh, pg_id, and block_id. ctrl_desc_field (str) = Field to propagate from ctrl wells to ctrl field in httr trt grp cmp document, defaults to "stype" (vector) = Additional modifiers to determine trt_grp_id, "bl" = append block ID, grp_id_opts "pg" = append plate group ID, "pl" = append plate ID, "vs" = use trt_name for trt_wells and ctrl_wells, applied in the order specified, defaults to simple use only trt_name from trt_wells (str) = Name of collection with individual well treatment data, default is "httr well" well (str) = Name of collection with treatment group data, default is "httr_trt_grp_cmp" trt_grp (bool) = If document with same trt_grp_id exists already, should it be replaced? rerun Default is False db insert (bool) = Whether to insert the new document into database at all, Default is True

Value

(list) = Document that was inserted into trt_grp collection

treatPairGroup

treatPairGroup Construct a document for httr_trt_grp_cmp collection corresponding to a specific pair of treatment names. Given a specific pair of trt_name values, plus optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This is initially designed for the pairs of BioSpyder QC reference samples now used in newer screens, but should work for other situations as well. Note, if filtering reduces either trt or ctrl group to < min_reps samples this treatment group will be skipped.

Description

treatPairGroup Construct a document for httr_trt_grp_cmp collection corresponding to a specific pair of treatment names.

Given a specific pair of trt_name values, plus optional filter criteria, construct an appropriate document and insert into httr_trt_grp_cmp. This is initially designed for the pairs of BioSpyder QC reference samples now used in newer screens, but should work for other situations as well. Note, if filtering reduces either trt or ctrl group to < min_reps samples this treatment group will be skipped.

treatPairGroup 63

Usage

```
treatPairGroup(
  DB = getOption("DB_NAME"),
  trt_name,
  ctrl_name,
  trt_type = "QC sample",
  ctrl_type = "QC sample",
  pg_id = "",
  media = "",
  timeh = 0,
  trt_src = "",
  ctrl_src = "",
  qc_flags = c("OK"),
  min_reps = 2,
  well = getOption("HTTR_WELL_NAME"),
  ...
)
```

Arguments

DB	name of database with httr_well and httr_trt_grp_cmp collections
trt_name	(string) = Match to trt_name in httr_well
ctrl_name	(string) = Match to trt_name in httr_well, use as control group
trt_type	(string) = Match to stype in httr_well, defaults to "QC sample"
ctrl_type	(string) = Match to stype in httr_well, defaults to "QC sample"
pg_id	(string) = Filter both trt and ctrl wells, match to pg_id in httr_well, defaults to None
media	(string) = Filter both trt and ctrl wells, match to media in httr_well, defaults to None
timeh	(string) = Filter both trt and ctrl wells, match to timeh in httr_well, defaults to None
trt_src	(string) = Filter trt wells, match to rna_src field
ctrl_src	(string) = Filter ctrl wells, match to rna_src field
qc_flags	(vector) = Filter both trt and ctrl wells, match to qc_flag field, defaults to OK only
min_reps	(int) = Minimum number of replicates in trt and ctrl groups, respectively - if either is less, will generate an output message and return empty dict
well	$(string) = Name of collection with individual well treatment data, default is "httr_well"$
	= All additional args passed to treatGroupFromSamples

Value

(list) = Document that was inserted into trt_grp collection

64 tss_deg_count

trtGrpWells

trtGrpWells - Extract the complete set of wells involved in a series of treatment groups

Description

trtGrpWells - Extract the complete set of wells involved in a series of treatment groups

Usage

```
trtGrpWells(trt_list, subset = "all")
```

Arguments

trt_list (*list of list*) = List of docs from httr_trt_grp_cmp (list) that are expected to all have the same control wells

subset (character) = Single character value indicating which subset of wells to extract.

"all" (default) retrieves both ctrl and trt wells, "ctrl" and "trt" extrct just those

well subsets respectively.

Value

(character vector) = Unique set of sample_ids for all or subset of wells in trt_list

tss_deg_count

 tss_deg_count : Function to simply count the number or fraction of probes with |l2fc| >= cutoff (with optional pvalue cutoff)

Description

tss_deg_count: Function to simply count the number or fraction of probes with ||2fc| >= cutoff (with optional pvalue cutoff)

```
tss_deg_count(
  degs,
  abs_l2fc,
  pval = NULL,
  l2fc_col = "log2FoldChange",
  pval_col = "pvalue",
  as_frac = FALSE
)
```

tss_12_norm 65

Arguments

degs	(<i>data.frame</i>) = Results from httr_deg for a single trt_grp_id and anl_name, required.
abs_12fc	(numeric) = Count probes with 12fc >= this cutoff, required.
pval	(numeric) = Additional threshold for p-value, ignore by default.
l2fc_col	(character) = Name of 12fc column, defaults to "log2FoldChange" to match DESeq2 output format
pval_col	(character) = Name of p-value column, defaults to "pvalue" to match DESeq2 output format
as_frac	(<i>logical</i>) = Whether or not to return result as a fraction of total probes in table, defaults to FALSE

Value

(numeric or integer) = Either the total number or fraction of probes passing the threshold criteria

```
tss\_12\_norm: \quad \textit{tss\_l2\_norm:} \quad \textit{Function to compute Euclidean (L2) norm for top/bottom N probes by l2fc}
```

Description

tss_12_norm: Function to compute Euclidean (L2) norm for top/bottom N probes by 12fc

Usage

```
tss_12_norm(
  degs,
  n,
  pval = NULL,
  12fc_col = "log2FoldChange",
  pval_col = "pvalue",
  rank_col = 12fc_col,
  rank_dir = "both",
  debug = getOption("debug", default = FALSE)
)
```

Arguments

degs	(<i>data.frame</i>) = Results from httr_deg for a single trt_grp_id and anl_name, required.
n	(integer) = Use the top and bottom n probes by 12fc
pval	(<i>numeric</i>) = If specified, change all l2fc to 0 for probes with p-values above this threshold before ranking, defaults to NULL
l2fc_col	(character) = Name of 12fc column, defaults to "log2FoldChange" to match DESeq2 output format
pval_col	(character) = Name of pvalue column, defaults to "pvalue" to match DESeq2 output format, only matters if pval is specified

66 unit_subtest_ident

rank_col	(character) = Which column to use for ranking the top/bottom N probes, defaults to 12fc_col but can rank by e.g. pval or SNR metric
rank_dir	$({\it character})$ = Whether to use the top N, bottom N, or both top/bottom N when computing L2 Norm, defaults to both
debug	(logical) = Whether to print debug messages

Value

(numeric) = The euclidean L2 norm of the top/bottom n l2fc values

```
unit_subtest_finish

unit_subtest_finish generate a single row of output for the outcome of

one test in a unit test and append it to current results
```

Description

unit_subtest_finish generate a single row of output for the outcome of one test in a unit test and append it to current results

Usage

```
unit_subtest_finish(col = NA, row = NA, test_name, test_status, input)
```

Arguments

```
col (int) = olumn number within the entire tested object that the test was run on
row (int) = column number within the entire tested object that the test was run on
test_name (character) = description of what was tested
test_status (character) = either "PASS" or "FAIL" - the outcome of the test
input (data.frame) = the growing result object report - 4 columns deep
```

Value

Value (data.frame) = returns a dataframe with 4 columns

```
unit_subtest_ident unit_subtest_set - test whether two vectors contain the same set of el-
ements
```

Description

unit_subtest_set - test whether two vectors contain the same set of elements

```
unit_subtest_ident(result_1, result_2, running_result, row = NA, col = NA)
```

unit_subtest_length 67

Arguments

```
result_1 (vector) = object to be tested

result_2 (vector) = object to be used as a reference

running_result
(data.frame) = result object to which subtest results should be appanded

row (int) = row number within the entire tested object that the test was run on

col (int) = column number within the entire tested object that the test was run on
```

Value

(data.frame) = returns a dataframe with 4 columns

```
unit_subtest_length

unit_subtest_length - test whether two objects are of the same length
```

Description

unit_subtest_length - test whether two objects are of the same length

Usage

```
unit_subtest_length(result_1, result_2, running_result, col = NA, row = NA)
```

Arguments

```
(data.frame) = returns a dataframe with 4 columns
```

68 unit_subtest_type

```
unit_subtest_set unit_subtest_set - test whether two vectors contain the same set of elements
```

Description

unit_subtest_set - test whether two vectors contain the same set of elements

Usage

```
unit_subtest_set(result_1, result_2, running_result, col = NA, row = NA)
```

Arguments

Description

 $unit_subtest_type \text{ - test } whether \text{ two objects are of the same type - a wrapper around } typeof()$

Usage

```
unit_subtest_type(result_1, result_2, running_result)
```

Arguments

```
(data.frame) = returns a dataframe with 4 columns
```

unit_subtest_values 69

```
unit\_subtest\_values \\ unit\_subtest\_values - depreciated - use \ unit\_subtest\_values\_1 \ instead
```

Description

```
unit_subtest_values - depreciated - use unit_subtest_values_1 instead
```

Usage

```
unit_subtest_values(result_1, result_2, running_result, tolerance)
```

```
unit\_subtest\_values\_1\\ unit\_subtest\_values\_1 - test\ whether\ two\ vectors\ contain\ the\ same\ values
```

Description

unit_subtest_values_1 - test whether two vectors contain the same values

Usage

```
unit_subtest_values_1(result_1, result_2, running_result, tolerance)
```

Arguments

```
(data.frame) = returns a dataframe with 4 columns
```

70 updateqcFlags

unit_test

unit_test- test whether two data frames contain the same values

Description

unit_test- test whether two data frames contain the same values

Usage

```
unit_test(result_1, result_2, tolerance)
```

Arguments

Value

(data.frame) = returns a dataframe with 4 columns

updateqcFlags

updateqcFlags Function to update QC_flag for samples based on changed QC_flag threshold. Updates are saved into httr_counts_qc collection function iterates through all documents found in httr_count_qc or iterates through all documents found in httr_count_qc and the result of a passed query to httr_well_trt (query passed in the ...)

Description

updateqcFlags Function to update QC_flag for samples based on changed QC_flag threshold. Updates are saved into httr_counts_qc collection function iterates through all documents found in httr_count_qc or iterates through all documents found in httr_count_qc and the result of a passed query to httr_well_trt (query passed in the ...)

```
updateqcFlags(
  db_host,
  db_name,
  db_collection = "httr_counts_qc",
  min_mapd_frac = get_global_option("min_mapd_frac"),
  min_n_reads_mapd = get_global_option("min_n_reads_mapd"),
  min_n_sig80 = get_global_option("min_n_sig80"),
  min_n_cov5 = get_global_option("min_n_cov5"),
  max_top10_prop = get_global_option("max_top10_prop"),
  max_gini_coef = get_global_option("max_gini_coef"),
```

```
qc_flags = get_global_option("qc_flags"),
  debug = getOption("debug", default = FALSE),
   ...
)
```

Arguments

```
db_host
                (character) = MongoDB host URL
                (character) = MongoDB name
db_name
db_collection
                 (character) = Name of counts_qc collection, default: httr_counts_qc
min mapd frac
                 (numeric) = minimum mapd_frac, flag with LOW_MAPD_FRAC below this
min_n_reads_mapd
                (integer) = minimum n_reads_mapd, flag with LOW_READS below this cutoff
                (integer) = minimum n_sig80, flag with LOW_NSIG80 below this cutoff
min_n_sig80
                (integer) = minimum n_cov5, flag with LOW_NCOV5 below this cutoff
min_n_cov5
max_top10_prop
                 (numeric) = maximum top10_prop, flag with HIGH_TOP10 above this cutoff
max_gini_coef
                (numeric) = maximum gini_coef, flag with HIGH_GINI above this cutoff
                (character vector) = list of flags to apply in priority order
qc_flags
                (logical) = Whether to report debug messages
debug
```

Value

nothing explicitly

```
update_attenuation_factors
```

update_attenuation_factors Create function in httr/Rlib/db/probe.R to update/replace attenuations factors in an existing httr_probe collection

Description

update_attenuation_factors Create function in httr/Rlib/db/probe.R to update/replace attenuations factors in an existing httr_probe collection

Usage

```
update_attenuation_factors(host, db, input_file, collection = "httr_probe")
```

Arguments

host: mongo url with or without port db: mongo database or sandbox

collection: mongo collection to update - default httr probe

changes attenuation values for probes found in dataframe in the collection and

change to '1' any probe in collection not in the data.frame

72 validateProbeManifest

Value

nothing explicit

validateProbeManifest

validateProbeManifest Function to validate a probe manifest (taken as a data.frame), so that there is a check to ensure the probe manifest is formatted correctly, is not missing and required data, and can be used.

Description

validateProbeManifest Function to validate a probe manifest (taken as a data.frame), so that there is a check to ensure the probe manifest is formatted correctly, is not missing and required data, and can be used.

Usage

```
validateProbeManifest(input_file)
```

Arguments

input_file = defaults to ../data/httrpl_automationTestData/httr_probe_data/human_wt_1.2.csv constraints are define @ https://teams.microsoft.com/l/entity/com.microsoft.teamspace.tab.wiki/tab::' f6ec-490c-87e1-685a2c988b22?context=%7B%22subEntityId%22%3A%22%7B%5C%22pageId% 6748-4867-acf9-76aacbeca6a7

Value

a list of validation errors or an empty list if all tests pass

Index

tests, 60 .onLoad, 4 getDBGs, 32 getDocIDs, 33 getFCmatrix, 34 getProbeManifest, 35 getTrts, 36 getWellCounts, 36 getWellCounts, 36 getWellCounts, 36 getWellCounts, 36 getWellInfo, 37 allRefRoups, 9 allTestGroups, 10 as.mongo.date, 11 bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countGcfini, 17 countQc, 18 countQcflag, 19 countStats, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt, 23 fill_dose_for_well_trt, 23 fill_dose_for_well_trt, 23 fill_dose_for_well_trt, 25 Find, 26 findByID, 26 findByID, 26 finddDs, 27 flattenDeG, 28 format_for_dbe, 16 get_global_option, 38 getFCmatrix, 34 getProbeManifest, 35 getDBGs, 32 getDocIDs, 33 getFCmatrix, 34 getProbeManifest, 35 getPcoletia, 34 getProbeManifest, 35 getVellouts, 36 getDBGs, 32 getDocIDs, 33 getFCmatrix, 34 getProbeManifest, 35 getCdefault, 35 getClefault, 35 getCdefault, 35 getCdefault, 35 getCdefault, 35 getCdefault, 35 getCdefault, 35 getClefault, 35 getCdefault, 35 getClefault, 35 getClefult, 36 getWellounts, 36 getWellounts, 36 getWellou	Tania datasats	
.onLoad, 4 addQCplotCols, 5 allBspLysateGroups, 6 allBspRNAGroups, 6 allBspRNAGroups, 7 allRefDoseGroups, 8 allRefGroups, 9 allRefGroups, 9 allTestGroups, 10 as.mongo.date, 11 bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 plotQCjitter, 47 plotQCviolin, 47 process_chem_columns, 49 process_chem_lolumns, 49 process_chem_lo	*Topic datasets	getCredentials, 31
addQCplotCols,5 allBspLysateGroups,6 allBspRNAGroups,6 allBspRNAGroups,7 allRefDoseGroups,7 allRefGroups,9 allRefRNAGroups,9 allRefRNAGroups,9 allRefRNAGroups,10 as.mongo.date,11 bspRNAGroup,12 bulkLysateGroup,14 check_collection_remapped,15 chemTreatGroup,16 countGini,17 countQC,18 countQCflag,19 countStats,20 degFrame, 20 deleteByID, 21 DESeq2_cut,22 extract_and_run,22 fill_dose_for_well_trt,23 fill_dose_for_well_trt, 23 fill_dose_for_well_trt, 23 fill_dose_for_well_trt, 25 Find, 26 findByID, 26 findByID, 26 findTbs, 27 flattenDEG, 28 format_for_dbs insert, 28 getDocIDs, 33 getFCmatrix, 34 getProbeManifest, 35 getPCmatrix, 34 getProbeManifest, 35 getQCdefault, 35 getWellCounts, 36 getDCcdefault, 35 getCdefault, 35 getColefault, 35 getColefault, 35 getColefault, 35 getClefault,		-
addCCplotCols,5 allBspRNaGroups,6 allBspRNaGroups,6 allBspRNaGroups,7 allRefDoseGroups,8 allRefGroups,9 allRefRnaGroups,9 allTestGroups,10 as.mongo.date,11 bspRNaGroup,12 bulkLysateGroup,14 check_collection_remapped,15 countGini,17 countQC,18 countQCflag,19 countStats,20 deleteByID,21 DESeq2_cut,22 fill_dose_for_well_trt_vapper,24 filterSamplesByTotal,25 Find,26 findByID,26 findByID,26 findIDs,27 flattenDEG,28 format_for_dbs. fill_dose_tor_dbs. format_for_dbs. format_for_dbs. format_for_dbs. format_for_dbs. getFCmatrix,34 getProbeManifest,35 getQCdefault,35 getQCdefault,35 getQCdefault,35 getQCdefault,35 getWellCounts,36 getWellCounts,36 getWellCounts,36 getWellCounts,36 getWellCounts,36 getWellInfo,37 allRefRnAGroup,59 insertLand_httr_study,41 insertByID,39 insertManyDEG,39 insertManyDEG,40 iterate,44 mongoQuery,44	.onLoad, 4	-
allBspLysateGroups, 6 allBspRNAGroups, 6 allBspRNAGroups, 6 allBukLysateGroups, 7 allRefDoseGroups, 8 allRefGroups, 9 allRefGroups, 9 allRefGroups, 10 as.mongo.date, 11 bspLysateGroup, 12 bulkLysateGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countQC, 18 countQC, 18 countQCflag, 19 countStats, 20 deleteByID, 21 bESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt_varapper, 24 filterProbesByMean, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_dby long getProbeManifest, 35 getQCdefault, 35 getWellTonts, 36 getWellCounts, 36 getWellCounts, 36 getWellCounts, 36 getWellCounts, 36 getWellCounts, 36 getWellInfo, 37 insert_into_httr_study, 41 insertByID, 39 insertConeDEG, 40 insert_oneDEG, 40 insert_oneDeG, 40 insert_monyDeG, 39 insertManyDEG, 39 insertManyDEG, 39 insertMonple, 44 mongoQuery, 44 mongoQuery, 44 mongoURL, 45 cuterFences, 46 outerFences, 56 outerFences, 46 outerFences, 56 outerFences, 46 outerFences, 49 ou	addOCnlotCols 5	-
allBspRNAGroups, 6 allBulkLysateGroups, 7 allRefDoseGroups, 8 allRefGroups, 9 allRefRnAGroups, 9 allTestGroups, 10 as.mongo.date, 11 bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt_wrapper, 24 filterSamplesByTotal, 25 Find, 26 findIDs, 27 filattenDEG, 28 format_for_dbyinsert, 28 getQCdefault, 35 getTrts, 36 getWellInfo, 37 getWellCounts, 36 getWellInfo, 37 getWellCounts, 36 getWellInfo, 37 getWellCounts, 36 getWellInfo, 37 getWellCounts, 36 getWellInfo, 37 insert_into_httr_study, 41 insertByID, 39 insertOneDEG, 40 iterate, 44 mongoUery, 44 mongoUery, 44 mongoURL, 45 counterFences, 46 counterfence, 49 counterfence, 49 counterfence, 40 counterfence, 49 counterfence, 40 coun		-
allBulkLysateGroups,7 allRefDoseGroups,8 allRefGroups,9 allRefRnAGroups,9 allTestGroups,10 as.mongo.date,11 bspLysateGroup,12 bulkLysateGroup,14 check_collection_remapped,15 chemTreatGroup,16 countGini,17 countQC,18 countQCflag,19 countStats,20 degFrame,20 deleteByID,21 DESeq2_cut,22 extract_and_run,22 fill_dose_for_well_trt,23 fill_dose_for_well_trt_wrapper, 24 filterProbesByMean, 24 filterSamplesByTotal,25 Find, 26 findByID, 26 findIDs,27 flattenDEG, 28 format_for_db_insert,28 getTrts, 36 getWellCounts, 36 getWellInfo, 37 allRefRoaGroups, 30 insert_into_httr_study, 41 insertByID, 39 insertManyDEG, 49 insertManyDEG, 39 insertManyDEG, 49 insertManyDEG, 40 insertManyDE, 39 insertManyDEG, 49 insertManyDE, 4 mongoult, 45 mongoult, 45 mongoult, 45 mongoult, 45 mongoult, 45 mongou		
allRefDoseGroups, 8 allRefGroups, 9 allRefRNAGroups, 9 allTestGroups, 10 as.mongo.date, 11 bspLysateGroup, 12 bulkLysateGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper, 24 filterProbesByMean, 24 filterSamplesByTotal, 25 Find, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 getWellCounts, 36 getWellInfo, 37 getWellCounts, 36 getWellCounts, 36 getWellInfo, 37 insert_into_httr_study, 41 insert_pite, 50 getWellInfo, 37 ansert_into_httr_study, 41 insert_pite, 50 getWellCounts, 36 getWellCounts, 36 getWellCounts, 36 getWellCounts, 36 getWellCounts, 36 getWellInfo, 37 insert_into_httr_study, 41 insert_pite, 41 insert_pite, 50 getWellInfo, 37 ansert_into_httr_study, 41 insert_pite, 51 insert_into_httr_study, 41 insert_pite, 51 onert_into_httr_study, 41 insert_pite, 52 onert_into_httr_study, 41 insert_pite, 52 onert_into_httr_study, 41 insert_pite, 52 onert_into_httr_study, 41 insert_pite, 52 onert_into_httr_study, 41 insert_pite, 47 onert_into_httr_study, 41 onert_into_httr_study, 41 onert_into_httr_study, 41 onert_into_httr_study, 41 onert_into_httr_study.		-
allRefGroups, 9 allRefRNAGroups, 9 allTestGroups, 10 as.mongo.date, 11 bspLysateGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt_wrapper, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findByID, 26 findByID, 26 findByID, 28 format_for_db_insert, 28 getWellInfo, 37 insertLinto_httr_study, 41 insertByID, 37 insertManyDEG, 39 insertOneDEG, 40 iterate, 44 bulkLysateGroup, 14 mongoQuery, 44 mongoURL, 45 copenMongo, 45 outerFences, 46 couterFences, 46 ropenMongo, 45 outerFences, 46 couterFences, 46 ropenMongo, 45 ropenM		-
allRefRNAGroups, 9 allTestGroups, 10 as.mongo.date, 11 bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt_wrapper, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findByID, 26 findByID, 26 findTestGroup, 18 countGoup, 29 get_global_option, 38 getAnlName, 29 insert_into_httr_study, 41 insertByID, 21 insertByID, 39 insertManyDEG, 39 insertOneDtEG, 40 insertManyDEG, 39 insertManyDEG, 40 iterate, 44 mongoURL, 45 couterFences, 46 couterFences, 47 mongoURL, 45 couter the MongoURL the MongoURL couter the MongoURL the Mongo		
allTestGroups, 10 as.mongo.date, 11 bspLysateGroup, 11 bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20		getwellinio, 3/
as.mongo.date, 11 bspLysateGroup, 11 bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper, 24 filterSamplesByMean, 24 filterSamplesByTotal, 25 Find, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 getAnlName, 29 insertByID, 39 insertManyDEG, 39 insertManyDEG, 39 insertManyDEG, 40 mongoQuery, 44 mongoURL, 45 countGille, 44 check_collectoup Fook and countGille, 45 countGille, 44 check_collectoup Fook and countGille, 45 countGille, 44 check_collectoup Fook and countGille, 45 countGille, 44 check_collectoup, 44 mongoURL, 45 countGille, 44 check_collectoup Fook and countGille, 45 countGille, 44 check_collectoup Fook and countGille, 45 countGille, 44 check_collectoup, 44 mongoURL, 45 countGille, 44 check_collectoup, 44 mongoURL, 45 countGille, 44 check_collectoup, 44 mongoURL, 45 countGille, 44 countGille, 44 countGille, 44 countGille, 44 countGille, 44 countGille, 45 countGille, 44 countGille, 44 countGille, 44 countGille, 44 countGille, 44 countGille, 45 countGille, 44 countGille, 45 countGille, 45 countGille, 45 coun	-	incort into http study 11
bspLysateGroup, 11 bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countQC, 18 countQC, 18 countQCflag, 19 countStats, 20 degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	_	
bspLysateGroup, 12 bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20	as.mongo.date, 11	——————————————————————————————————————
bspRNAGroup, 12 bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper, 24 filterProbesByMean, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 getAnlName, 29 iterate, 44 mongoQuery, 44 mongoURL, 45 countQCflag, 19 countStats, 40 plotQCjitter, 47 plotQCviolin, 47 process_chem_columns, 49 process_chem_Info, 48 process_chemInfo, 48 readWellTrtFile, 50 refChemGroup, 51 refRNAGroup, 52 runDESeq2, 53 runDESeq2, 53 runDESeq2, 53 runDESeq2Single, 56 sameCtrlWells, 57 sampleID_scan_and_update, 58 sampleID_wrapper, 59 splitIDs, 60 treatGroupFromSamples, 61	bspLysateGroup, 11	-
bulkLysateGroup, 14 check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 filterProbesByMean, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_dbe_insert, 28 getAnlName, 29 mongoQuery, 44 mongoURL, 45 countGini, 47 copenMongo, 45 couterFences, 46 couterFences, 4		
check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt_wrapper, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 getAnlName, 29 mongoQuery, 44 mongoURL, 45 openMongo, 45 outerFences, 46 ou		iterate, 44
check_collection_remapped, 15 chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper, 24 filterProbesByMean, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 getAnlName, 29 mongoURL, 45 mongoUnt, 46 mongoUnt, 45 mongoUnt, 45 mongoUnt, 45 mongoUnt, 45 mongoUnt, 40 mongoUnt, 45 mongoUnt, 46 mongoUnt, 47 mongoUnt, 47 mongoUnt, 47 mongoUnt, 47 m	, , , , , , , , , , , , , , , , , , ,	mongoOuery 44
chemTreatGroup, 16 countGini, 17 countQC, 18 countQCflag, 19 countStats, 20 plotQCjitter, 47 plotQCviolin, 47 process_chem_columns, 49 process_chemInfo, 48 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	check_collection_remapped, 15	
countQC, 18 countQCflag, 19 countStats, 20 plotQCjitter, 47 plotQCviolin, 47 process_chem_columns, 49 process_chemInfo, 48 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,		mongoons, 10
countQC, 18 countQCflag, 19 countStats, 20 plotQCjitter, 47 plotQCviolin, 47 degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	countGini, 17	openMongo, 45
countQCflag, 19 countStats, 20 plotQCyiolin, 47 plotQCviolin, 47 process_chem_columns, 49 process_chemInfo, 48 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	countQC, 18	
degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	countQCflag, 19	,
degFrame, 20 deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	countStats, 20	plotQCjitter,47
deleteByID, 21 DESeq2_cut, 22 extract_and_run, 22 extract_and_run, 22 readWellTrtFile, 50 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,		plotQCviolin,47
DESeq2_cut, 22 extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,		process_chem_columns,49
qcBatch, 49 extract_and_run, 22 readWellTrtFile, 50 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	-	process_chemInfo,48
extract_and_run, 22 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	DESeq2_cut, 22	
readWellTrtFile, 50 fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,		qcBatch, 49
fill_dose_for_well_trt, 23 fill_dose_for_well_trt_wrapper,	extract_and_run, 22	
fill_dose_for_well_trt_wrapper, 24 filterProbesByMean, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 get_global_option, 38 getAnlName, 29 refRNAGroup, 52 runDESeq2, 53 runDESeq2ForChemCond, 54 runDESeq2Single, 56 sameCtrlWells, 57 sampleID_scan_and_update, 58 sampleID_wrapper, 59 splitIDs, 60 tests, 60 treatGroupFromSamples, 61	fill dose for well tot 23	
runDESeq2,53 filterProbesByMean, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 get_global_option, 38 getAnlName, 29 runDESeq2, 53 runDESeq2ForChemCond, 54 run		-
filterProbesByMean, 24 filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 get_global_option, 38 getAnlName, 29 runDESeq2ForChemCond, 54 runDESeq2Single, 56 sameCtrlWells, 57 sampleID_scan_and_update, 58 sampleID_wrapper, 59 splitIDs, 60 tests, 60 treatGroupFromSamples, 61		
filterSamplesByTotal, 25 Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 get_global_option, 38 getAnlName, 29 runDESeq2Single, 56 sameCtrlWells, 57 sampleID_scan_and_update, 58 sampleID_wrapper, 59 splitIDs, 60 tests, 60 treatGroupFromSamples, 61		=
Find, 26 findByID, 26 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 get_global_option, 38 getAnlName, 29 sameCtrlWells, 57 sampleID_scan_and_update, 58 sampleID_wrapper, 59 splitIDs, 60 tests, 60 treatGroupFromSamples, 61		-
findByID, 26 findIDs, 27 findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 get_global_option, 38 getAnlName, 29 sameCtrlWells, 57 sampleID_scan_and_update, 58 sampleID_wrapper, 59 splitIDs, 60 tests, 60 treatGroupFromSamples, 61		runDESeq2Single, 36
findIDs, 27 flattenDEG, 28 format_for_db_insert, 28 get_global_option, 38 getAnlName, 29 sampleID_scan_and_update, 58 sampleID_wrapper, 59 splitIDs, 60 tests, 60 treatGroupFromSamples, 61		sameCtrlWells 57
flattenDEG, 28 sampleID_wrapper, 59 splitIDs, 60 set_global_option, 38 tests, 60 getAnlName, 29 treatGroupFromSamples, 61		
format_for_db_insert, 28 splitIDs, 60 get_global_option, 38 tests, 60 getAnlName, 29 treatGroupFromSamples, 61		
get_global_option, 38 tests, 60 getAnlName, 29 treatGroupFromSamples, 61		
getAnlName, 29 treatGroupFromSamples, 61	TOTHLAC_TOT_UD_THISETC, 20	3P11C1D3, 00
getAnlName, 29 treatGroupFromSamples, 61	get_global_option, 38	tests, 60
		treatPairGroup, 62

74 INDEX

```
trtGrpWells, 64

tss_deg_count, 64

tss_12_norm, 65

unit_subtest_finish, 66
unit_subtest_ident, 66
unit_subtest_length, 67
unit_subtest_set, 68
unit_subtest_type, 68
unit_subtest_values, 69
unit_subtest_values_1, 69
unit_test, 70
update_attenuation_factors, 71
updateqcFlags, 70

validateProbeManifest, 72
```